


DEFENDER[™] SERIES SCALE ASSEMBLY INSTRUCTIONS



P/N 80252786 © 2009 Ohaus Corporation, all rights reserved.

Printed in China





3000 Series Indicators Instruction Manual





T31XW Indicator

TABLE OF CONTENTS

 1.1 Safety Precautions	EN-5 EN-9
1.3 Control Functions	EN-9
	N-10
2. INSTALLATIONE	11-10
2.1 UnpackingE	N-10
2.2 External ConnectionsE	
2.2.1 RS232 Interface Cable to T31PE	N-10
2.2.2 AC Power to T31PE	N-10
2.2.3 AC Power to T31XWE	N-10
2.2.4 Battery power (T31P)E	N-11
2.2.5 Mounting Bracket to T31XWE	N-11
2.3 Internal ConnectionsE	N-11
2.3.1 Opening the HousingE	
2.3.2 Scale Base to T31P or T31XW	
2.3.3 RS232 Interface Cable to T31XWE	N-12
2.4 T31P Rear Cover Orientation	N-13
2.5 Direct Wall Mounting (T31P only)E	
2.6 Mounting Bracket (T31XW only)E	N-13
3. SETTINGSE	N-14
3.1 Menu Structure	N-14
3.2 Menu NavigationE	N-15
3.3 Calibration MenuE	N-15
3.3.1 Span CalibrationE	N-16
3.3.2 Linearity CalibrationE	N-16
3.3.3 Geographical Adjustment FactorE	N-17
3.3.4 End CalibrationE	N-17
3.4 Setup MenuE	N-19
3.4.1 ResetE	N-19
3.4.2 Legal for tradeE	
3.4.3 Calibration UnitE	
3.4.4 CapacityE	
3.4.5 GraduationE	
3.4.6 Power On UnitE	
3.4.7 Zero RangeE	
3.4.8 End Setup	

TABLE OF CONTENTS (Cont.)

3.5	Readou	ıt Menu	.EN-21
	3.5.1	Reset	.EN-21
	3.5.2	Filter	.EN-22
	3.5.3	Auto-Zero Tracking	.EN-22
	3.5.4	Backlight	.EN-22
	3.5.5	Auto Off Timer	.EN-22
	3.5.6	End Readout	.EN-22
3.6		Menu	
	3.6.1	Reset	.EN-23
	3.6.2	Parts Counting Mode	.EN-23
	3.6.3	End Mode	.EN-23
3.7	Unit Me	ะกน	.EN-24
	3.7.1	Reset	.EN-24
	3.7.2	Kilogram Unit	.EN-24
	3.7.3	Pound Unit	.EN-24
	3.7.4	Gram Unit	.EN-24
	3.7.5	Ounce Unit	.EN-24
	3.7.6	Pound Ounce Unit	.EN-24
	3.7.7	End Unit	.EN-24
3.8	Print M	enu	.EN-25
		Reset	
		Baud	
	3.8.3	Parity	
	3.8.4	Stop Bit	
		Handshake	
		Print Stable Data Only	
		Auto Print	
		Content	
		End Print	
3.9		ock Menu	
	3.9.1	Reset	.EN-27
	3.9.2	Lock Calibration	.EN-27
	3.9.3	Lock Setup	.EN-27
	3.9.4	Lock Readout	.EN-27
	3.9.5	Lock Mode	.EN-27
	3.9.6	Lock Unit	.EN-27
	3.9.7	Lock Print	.EN-28
	3.9.8	End Lock	.EN-28

TABLE OF CONTENTS (Cont.)

3.10	O Security Switch	EN-28
4.	OPERATION	EN-28
4.1	Turning Indicator On/Off	EN-28
4.2	Zero Operation	EN-28
4.3	Manual Tare	EN-28
4.4	Changing Units of Measure	EN-29
4.5	Printing Data	EN-29
4.6	Application Modes	EN-29
	4.6.1 Weighing	EN-29
	4.6.2 Parts Counting	EN-29
5.	SERIAL COMMUNICATION	
	Interface Commands	
	Output Format	
6.	LEGAL FOR TRADE	EN-33
	Settings	
	Verification	
	Sealing	
0.0	6.3.1 Physical Seals	
	6.3.2 Audit Trail Seal	
7.	MAINTENANCE	EN-36
7.1	Model T31P Cleaning	
7.2	Model T31XW Cleaning	EN-36
7.3	Troubleshooting	EN-34
7.4	Service Information	EN-37
8.	TECHNICAL DATA	EN-38
8.1	Specifications	EN-38
8.2	Accessories and Options	EN-39
8.3	Drawings and Dimensions	EN-40
8.4	Compliance	EN-41

1. INTRODUCTION

This manual contains installation, operation and maintenance instructions for the T31P and T31XW Indicators. Please read this manual completely before installation and operation.


1.1 Safety Precautions



For safe and dependable operation of this equipment, please comply with the following safety precautions:

- Verify that the input voltage range printed on the data label matches the local AC power to be used.
- Make sure that the power cord does not pose a potential obstacle or tripping hazard.
- Use only approved accessories and peripherals.
- Operate the equipment only under ambient conditions specified in these instructions.
- Disconnect the equipment from the power supply before cleaning.
- Do not operate the equipment in hazardous or unstable environments.
- Do not immerse the equipment in water or other liquids.
- Service should only be performed by authorized personnel.
- The T31XW is supplied with a grounded power cable. Use only with a compatible grounded power outlet.

1.2 Overview of Parts and Controls







Figure 1-1. T31P Indicator.

TABLE	1-1.	T31P	PARTS.
-------	------	-------------	--------


ltem	Description	
1	Data Label	
2	Front Housing	
3	Control Panel	
4	Security Screw	
5	Key Hole (4) for wall	
	mounting	
6	Screw (4)	
7	Data Label	
8	Rear Housing	
9	Power Receptacle	
10	Strain Relief for Load Cell	
	Cable	
11	RS232 Connector	



1.2 Overview of Parts and Controls (Cont.)

TABLE 1-2. T31XW PARTS.

Item	Description	
1	Data Label	
2	Front Housing	
3	Adjusting Knob (2)	
4	Control Panel	
5	Mounting Bracket	
6	Screw (4)	
7	Rear Housing	
8	Data Label	
9	Strain Relief for RS232	
10	Strain Relief for Load Cell	
	Cable	
11	Power cord	







Figure 1-2. T31XW Indicator.

1.2 Overview of Parts and Controls (Cont.)





LOAD CELL WIRING





TABLE 1-3. MAIN PC BOARD.

Item	Description	
1	Keypad Connector J4 T31XW Model only	
2	Battery Connector (T31P only)	
3	Line Power Input	
4	Sense Jumper W1	
5	Load Cell Terminal Block J5	
6	Sense Jumper W2	
7	RS232 Terminal Block J7 T31XW Model only	
8	RS232 Connector J6 T31P Model only	
9	LFT On / Off Switch	

RS232 WIRING



1.2 Overview of Parts and Controls (Cont.)




Figure 1-4. Controls and Indicators.

No.	Designation	
1	Capacity Label Window	
2	TARE <i>Menu</i> button	
3	TARE symbol	
4	FUNCTION Mode button	
5	NET symbol	
6	PCS symbol	
7	PRINT Units button	
8	Calibration Mode symbol	
9	Battery symbol (T31P only)	
10	ON/ZERO Off button	
11	Center of Zero symbol	

TABLE 1-4. CONTROL PANEL.

1.3 Control Functions

Button	On/Zero Off Yes	Print Units No	Function Mode Back	Tare Menu Exit
Primary Function	ON/ZERO	PRINT	FUNCTION	TARE
(Short Press)	If Indicator is On, sets zero.	Sends the current value to the COM port if AUTOPRINT is set to Off.	Initiates an application mode.	Performs a tare operation.
Secondary Function	Off	Units	Mode	Menu
	-			
(Long Press)	Turns the Indicator on or off.	Changes the weighing Unit.	Allows changing the application mode.	Enter the User menu.
				View the Audit Trail event
			Press and hold allows	counters (extended
			scrolling through modes.	press)
Menu Function	Yes	No	Back	Exit
(Short Press)	Accepts the current	Advances to the next	Moves Back to previous	Exits the User menu.
	setting on the display.	menu or menu item.	menu item.	
				Aborts the calibration in
		Rejects the current	Decrements the value.	progress.
		setting on the display		
		and advances to the		
		next available setting.		
		Increments the value.		

TABLE 1-5. CONTROL FUNCTIONS.

2. INSTALLATION

2.1 Unpacking

Unpack the following items:

- T31P or T31XW Indicator
- AC Adapter (T31P only)
- Mounting Bracket (supplied with T31XW only)
- Knobs (2) (supplied with T31XW only)
- Capacity Label Sheet
- Instruction Manual CD
- Warranty Card
- LFT sealing Kit

2.2 External Connections

2.2.1 RS232 interface Cable to T31P

Connect the optional RS232 cable to the RS232 connector Figure 1-1, item 13).

Pin	Connection
1	N/C
2	TXD
3	RXD
4	N/C
5	GND
6	N/C
7	N/C
8	N/C
9	N/C




Figure 2-1. RS232 Pins.

2.2.2 AC Power to T31P

Connect the AC Adapter to the power receptacle (Figure 1-1, item 8), then plug the AC Adapter into an electrical outlet.

2.2.3 AC Power to T31XW

Connect the AC plug to a properly grounded electrical outlet.

2.2.4 Battery Power (T31P Only)

The indicator can be operated on the internal rechargeable battery when AC power is not available. The indicator will automatically switch to battery operation if there is a power failure or the power cord is removed.


Note:

Before using the indicator for the first time, the internal rechargeable battery should be fully charged for up to 12 hours. The indicator can be operated during the charging process. The battery is protected against over charging and the indicator can remain connected to the AC power line.

Connect AC power to the indicator and allow it to charge. While the battery is charging, the triangle above the battery function symbol will light. When the battery is fully charged, this triangle will disappear.

The indicator can operate for up to 100 hours on a fully charged battery.

During battery operation, a flashing triangle above the battery function symbol indicates the battery is low and requires recharging. Approximately 60 minutes of operation will remain when the battery symbol starts to blink. The indicator will display Lo.BAT and automatically turn off when the battery is fully discharged.



2.2.5 Mounting Bracket to T31XW

Align the mounting bracket over the threaded holes in the side of the indicator and install the knobs. Adjust the indicator to the desired angle and tighten the knobs.

2.3 Internal Connections

Some connections require the housing to be opened.

2.3.1 Opening the Housing



CAUTION: ELECTRICAL SHOCK HAZARD. REMOVE ALL POWER CONNECTIONS TO THE INDICATOR BEFORE SERVICING OR MAKING INTERNAL CONNECTIONS. THE HOUSING SHOULD ONLY BE OPENED BY AUTHORIZED AND QUALIFIED PERSONNEL, SUCH AS AN ELECTRICAL TECHNICIAN.

EN-12

T31P

Remove the four Phillips head screws from the rear housing. Open the housing being careful not to disturb the internal connections. Once all connections are made, reattach the front housing.

T31XW

Remove the four hex head screws from the rear housing. Open the housing by carefully pulling the top of the front housing forward. Once all connections are made, reattach the front housing. The screws should be tightened fully to maintain a watertight seal.

2.3.2 Scale Base to T31P or T31XW


Pass the load cell cable through the strain relief (Figure 1-1, item 9 or Figure 1-2, item 10) and attach it to terminal block J5 (Figure 1-3, item 5). Re-tighten the strain relief to ensure a watertight seal.

Jumper Connections

For a 4-wire load cell with no sense wires: Jumpers W2 and W3 must be shorted.

For a 6-wire load cell that includes sense wires, see Figure 2-2. Jumpers W2 and W3 must be opened.

For load cells with an extra ground shield wire: Connect the shield to the center position (GND) of J5.





After wiring is completed and jumpers are in place, replace the indicator housing screws. Make sure the strain relief is properly tightened.

2.3.3 RS232 Interface Cable to T31XW

Pass the optional RS232 cable through the strain relief (Figure 1-2, item 9) and attach it to terminal block J7 (Figure 1-3, item 7). Re-tighten the strain relief to ensure a water tight seal.

Pin	Connection
J7-1	TXD
J7-2	RXD
J7-3	GND

Pin	Connection
J5-1	+EXCITATION
J5-2	+SENSE
J5-3	+SIGNAL
J5-4	GND
J5-5	-SIGNAL
J5-6	-SENSE
J5-7	-EXCITATION

2.4 T31P Rear Cover Orientation

The T31P is delivered in the wall mount orientation with the connections exiting below the display. The rear housing may be reversed so the connections exit above the display when the T31P is placed horizontally on a bench. See Figure 2-4. To reverse the rear housing, remove the four Phillips head screws, carefully rotate the housing 180°, and reinstall the screws.

CAUTION: Take care not to pinch any internal cables attached inside.


Figure 2-3. Wall Mount Configuration.



2.5 Direct Wall Mounting (T31P only)

The T31P indicator may be mounted directly to a wall using two screws (not included). Select appropriate size screws that fit into the holes at the bottom of the indicator housing. See Figure 2-5. When mounting to a wall without a solid backing, use appropriate anchoring hardware.

120 mm 4.7 in



30 mm

Figure 2-5. T31P Direct Wall Mounting.

2.6 Mounting Bracket (T31XW only)

Attach the bracket to a wall or table using fasteners (not supplied) that are appropriate for the type of mounting surface. The bracket will accommodate up to 6 mm (1/4") diameter screws. Locate the mounting holes as shown in Figure 2-6.





Figure 2-6. Mounting Bracket Dimensions.



EN-14

3 SETTINGS

3.1 Menu Structure

CALIBRATION-	> SETUP -	READOUT -	• MODE -	→ UNIT	→ PRINT	\rightarrow MENU LOCK \rightarrow END
→ SPAN	⊢RESET	→ RESET	→ RESET	-→ RESET	→ RESET	→ RESET
\rightarrow LINEARITY	⊢NO	Ь NO	Ь NO	⊢NO	⊢NО	⊢NO
-→ GEO	-→ YES	→ YES	→ YES	-→ YES	-→ YES	-→ YES
→ 0…31	→ LEGAL FOR TRADE	→ AVERAGING	-→ COUNT	→ KILOGRAM	⊢BAUD	⊢ LOCK CAL
\hookrightarrow END CAL	└→ OFF	LOW	└→ OFF	└→ OFF	→ 300, …19200	→ OFF
	Ь ON	→ MEDIUM	Ь ON	Ь ON	→ PARITY	Ч ON
	→ CALIBRATION UNIT	Ь НI	\hookrightarrow END MODE	→ POUND	→ 7 EVEN	\hookrightarrow LOCK SETUP
	-→ KILOGRAM	→ AUTO ZERO		-→ OFF	→7 ODD	GFF
	-→ POUND	└→ OFF		⊢ ОN	→ 7 NONE	⊔ ON
	→ CAPACITY	→0.5d		GRAM	→ 8 NONE	→ LOCK READOUT
	→ 5…20000	⊔ 1d		-→ OFF	└→ STOP	GFF
	→ GRADUATION	→ 3d		⊢ ОN	→ 1	Ч ON
	→ 0.001…20	\hookrightarrow EXPAND MODE		-→ OUNCE	→ 2	\hookrightarrow LOCK MODE
	→ POWER ON UNIT	-→ OFF		-→ OFF	HANDSHAKE	GFF
	-→ AUTO	⊢ ON		⊢ ОN	→ OFF	⊔ ON
	GRAM	→ BACKLIGHT		\mapsto POUND OUNCE	→ XON-XOFF	⊢ LOCK UNIT
	⊢KILOGRAM	→ AUTO		-→ OFF	STABLE ONLY	GFF
	-→ POUND	⊢ ON		⊢ ОN	⊢ OFF	Ь ON
	→OUNCE	└→ OFF		→ END UNIT	⊢ ON	→ LOCK PRINT
	\hookrightarrow POUND OUNCE	⊢AUTO OFF			→ AUTO PRINT	GFF
	→ ZERO RANGE	└→ OFF			Ь OFF	⊔ ON
	└→ 0%	⊢ SET 1			→ WHEN STABLE	→ END MENU LOCK
	₩2%	⊢ SET 2			→ INTERVAL	
	└→ 100%	⊢ SET 5			→ 1…3600	
	\mapsto END SETUP	\hookrightarrow END READOUT			→ CONTINUOUS	
					└→ CONTENT	
					⊢ GROSS	
					-→ NET	
					L→ TARE	
					-→ UNIT	
					\hookrightarrow END PRINT	

TABLE 3-1. MENU STRUCTURE.

3.2 Menu Navigation

TO ENTER THE MENU MODE

Press and hold the Menu button until MENU appears on the display. The first upper level menu appears on the display. Summary of button navigation functions in menu mode:

- --Yes Allows entry into the displayed menu.
 - Accepts the displayed setting and advances to the next menu item.
- --No Skips by the displayed menu.
 - Rejects the displayed setting or menu item and advances to the next available item.
- --Back Moves backwards through the upper and middle level menus.
 - Backs out of a list of selectable items to the previous middle level menu.
- --Exit Exits from menu directly to the active weighing mode.

3.3 Calibration Menu

Two calibration processes are available: Span Calibration and Linearity Calibration.

NOTES:

- 1. Make sure that appropriate calibration masses are available before beginning calibration.
- 2. Make sure that the scale base is level and stable during the entire calibration process.
- 3. Calibration is unavailable with LFT set to On.
- 4. Allow the Indicator to warm up for approximately 5 minutes after stabilizing to room temperature.
- 5. To abort calibration, press the **Exit** button anytime during the calibration process.

Span	Perform
Linearity	Perform
Geographic	
Adjustment	Set 00Set 19 Set 31
End Calibration	Exit CALIBRATE menu

Span Calibration uses two points to adjust the scale. The first point is the zero value where there is no weight on the scale. The second point is the Span value where a calibration mass is placed on the scale.

When SPAN is displayed, press the Yes button to access the Span Calibration menu item.

The display flashes 0.

With no weight on the scale, press the Yes button to establish the zero point.

The display shows --C-- while the zero point is established.

The display flashes the span calibration point. Place the specified weight on the scale and press the **Yes** button.

To choose a different span point, repeatedly press the **No** button to increment the selections or press the **Back** button to decrement the selections. Refer to Table 3-3 for available span points. When the desired value is displayed, place the specified weight on the scale and press the **Yes** button.

The display shows --C-- while the span point is established.

If span calibration was successful, the scale exits to the active weighing mode and displays the actual weight value.

3.3.2 Linearity Calibration

Linearity calibration uses 3 calibration points. The first calibration point is established with no weight on the scale. The second calibration point is established at approximately half capacity. The third calibration point is established at capacity. The Linearity calibration points are fixed and cannot be altered by the user during the calibration procedure. Refer to Table 3-3 for the linearity points.

When LINEAr is displayed, press the Yes button to access the Linearity Calibration menu item.

The display flashes 0. With no weight on the scale, press the Yes button to establish the zero point.

The display shows --C-- while the zero point is established.

The display flashes the mid calibration point.

Place the specified weight on the scale and press the Yes button.

The display shows --C-- while the mid point is established.

The display flashes the full calibration point.

Place the specified weight on the scale and press the Yes button.

The display shows --C-- while the full point is established.

If linearity calibration was successful, the scale exits to the active weighing mode and displays the actual weight value.




SPRA



- - [- -



∂ S kg		
ן באו		
	-	ןינס



L INERr

₩g kg
[

	∫ ∫ kg
	[
	J[] kg
	[
•	30.000 kg

3.3.3 Geographical Adjustment Factor

The Geographcial Adjustment Factor (GEO) is used to compensate for variations in gravity.

Note: Changing the GEO Factor alters the calibration. The GEO value was set at the factory and should only be changed by an authorized manufacturer's representative or certified verirication personnel.

Refer to table 3-2 to determine the GEO factor that corresponds to your location.

3.3.4 End Calibration

Advance to the next menu.

680

End

		Elevation in meters										
		0	325	650	975	1300	1625	1950	2275	2600	2925	3250
		325	650	975	1300	1625	1950	2275	2600	2925	3250	3575
						Ele	vation in t	feet				
		0	1060	2130	3200	4260	5330	6400	7460	8530	9600	10660
		1060	2130	3200	4260	5330	6400	7460	8530	9600	10660	11730
	tude						GEO value					
0°00′	5°46′	5	4	4	3	3	2	2	1	1	0	0
5°46′	9°52′	5	5	4	4	3	3	2	2	1	1	0
9°52′	12°44′	6	5	5	4	4	3	3	2	2	1	1
12°44′	15°06′	6	6	5	5	4	4	3	3	2	2	
15°06′	17°10′	7	6	6	5	5	4	4	3	3	2	2
17°10′	19°02′	7	7	6	6	5	5	4	4	3	3	2
19°02′	20°45′	8	7	7	6	6	5	5	4	4	3	3
20°45′	22°22′	8	8	7	7	6	6	5	5	4	4	3
22°22′	23°54′	9	8	8	7	7	6	6	5	5	4	4
23°54′	25°21′	9	9	8	8	7	7	6	6	5	5	4
25°21′ 26°45′	26°45′ 28°06′	10 10	9 10	9 9	8 9	8 8	7 8	7	6 7	6 6	5 6	5 5
26°45 28°06′	28°06 29°25′	10	10	10	9	9	8	8	7	6 7	6	5 6
28°06 29°25′	<u>29°25</u> 30°41′	11	10	10	10	9	9	8	8	7	6 7	6
<u>29 25</u> 30°41′	30 41 31°56′	11	11	10	10	10	9	9	8	8	7	7
31°56′	33°09′	12	12	11	10	10	10	9	9	8	8	7
33°09′	34°21′	12	12	12	11	10	10	10	9	9	8	8
34°21′	35°31′	13	13	12	12	11	11	10	10	9	9	8
35°31′	36°41′	14	13	13	12	12	11	11	10	10	9	9
36°41′	37°50′	14	14	13	13	12	12	11	11	10	10	9
37°50′	38°58′	15	14	14	13	13	12	12	11	11	10	10
38°58′	40°05′	15	15	14	14	13	13	12	12	11	11	10
40°05′	41°12′	16	15	15	14	14	13	13	12	12	11	11
41°12′	42°19′	16	16	15	15	14	14	13	13	12	12	11
42°19′	43°26′	17	16	16	15	15	14	14	13	13	12	12
43°26′	44°32′	17	17	16	16	15	15	14	14	13	13	12
44°32′	45°38′	18	17	17	16	16	15	15	14	14	13	13
45°38′	46°45′	18	18	17	17	16	16	15	15	14	14	13
46°45′	47°51′	19	18	18	17	17	16	16	15	15	14	14
47°51′	48°58′	19	19	18	18	17	17	16	16	15	15	14
48°58′	50°06′	20	19	19	18	18	17	17	16	16	15	15
50°06′	51°13′	20	20	19	19	18	18	17	17	16	16	15
51°13′	52°22′	21	20	20	19	19	18	18	17	17	16	16
52°22′	53°31′	21	21	20	20	19	19	18	18	17	17	16
53°31′	54°41′	22	21	21	20	20	19	19	18	18	17	17
54°41′	55°52′	22	22	21	21	20	20	19	19	18	18	17
55°52′	57°04′	23	22	22	21	21	20	20	19	19	18	18
57°04′	58°17′	23	23	22	22	21	21	20	20	19	19	18
58°17′	59°32′	24	23	23	22	22	21	21	20	20	19	19
59°32′	60°49′	24	24	23	23	22	22	21	21	20	20	19
60°49′	62°90′	25 25	24	24	23	23	22 23	22	21 22	21 21	20	20
62°90′ 63°30′	63°30′ 64°55′	25	25 25	24 25	24 24	23 24	23	22 23	22	21	21 21	20 21
63°55′	66°24′	26	25	25 25	24 25	24	23	23	22	22	21	21
66°24′	67°57′	20	26	25	25	24	24	23	23	22	22	21
67°57′	69°35′	27	20	26	26	25	24	24	23	23	22	22
69°35′	71°21′	28	27	20	26	26	25	24	24	23	23	22
71°21′	73°16′	28	27	27	20	26	26	25	24	24	23	23
73°16′	75°24′	20	28	28	27	20	26	26	25	25	24	23
75°24′	77°52′	29	20	28	28	27	20	26	26	25	25	24
77°52′	80°56′	30	29	29	28	28	27	27	26	26	25	25
80°56′	85°45′	30	30	29	29	28	28	27	27	26	26	25
85°45′	90°00′	31	30	30	29	29	28	28	27	27	26	26
00 -10	00 00			00	20	20	20	20	~ /	~ ~ /	20	20

TABLE 3-2. GEOGRAPHICAL ADJUSTMENT VALUES

3.4 Setup Menu

When the Indicator is used for the first time, enter this menu to set the Capacity and Graduation.

Reset	No, Yes
Legal for Trade	Off , On
Cal Unit	kg , lb
Capacity	520000
Graduation	0.001 20
Power On Unit	g, kg, lb, oz, lb:oz, Auto
Zero Range	0%, 2% , 100%
End Setup	Exit SETUP menu

3.4.1 Reset

Reset the Setup menu to the factory defaults.

No = not reset. Yes = reset.

NOTE: If the Legal for Trade menu item is set to ON, the Capacity, Graduation, Zero Range and Legal For Trade settings are not reset.

3.4.2 Legal for Trade

Set the legal for trade status.

OFF = off ON = on

Turning on the "LFT" menu setting has the following effects:

- Zero-range is set and locked on "2".
- Auto Zero Tracking is set and locked on 0.5d
- The lb:oz unit is not available as a power-on setting.

3.4.3 Calibration Unit

Set the unit during calibration.

- CAL UN kg = Calibrate using kg weights
- CAL UN Ib = Calibrate using pound weights

3.4.4 Capacity

Set the scale capacity from 5 to 20000. Refer to the Setup Table 3.3 for available settings.

58ευΡ

rESEE
00
985

LFE
OFF
00

ERL	.UN



TABLE 3-3. SETUP AND CALIBRATION VALUES

Capacity	ty Graduation size Graduation size with LFT ON with LFT OFF		Span calibration points	Linearity calibration points	
5	0.0005, 0.001, 0.002, 0.001, 0.002, 0.005		5	2, 5	
10	0.0005, 0.001, 0.002, 0.005, 0.01	0.002, 0.005, 0.01	5, 10	5, 10	
15 20	0.001, 0.002, 0.005, 0.01 0.001, 0.002, 0.005, 0.01, 0.02	0.005, 0.01 0.005, 0.01, 0.02	5, 10, 15 5, 10, 15, 20	5, 15 10, 20	
25	0.002, 0.005, 0.01, 0.02	0.005, 0.01, 0.02	5, 10, 15, 20, 25	10, 25	
30	0.002, 0.005, 0.01, 0.02	0.005, 0.01, 0.02	5, 10, 15, 20, 25, 30	15, 30	
40	0.002, 0.005, 0.01, 0.02	0.01, 0.02	5, 10, 15, 20, 25, 30, 40	20, 40	
50	0.005, 0.01, 0.02, 0.05	0.01, 0.02, 0.05	5, 10, 15, 20, 25, 30, 40, 50	25, 50	
60	0.005, 0.01, 0.02, 0.05	0.01, 0.02, 0.05	5, 10, 15, 20, 25, 30, 40, 50, 60	30, 60	
75	0.005, 0.01, 0.02, 0.05	0.02, 0.05	5, 10, 15, 20, 25, 30, 40, 50, 60, 75	30, 75	
100	0.005, 0.01, 0.02, 0.05,	0.02, 0.05, 0.1	5, 10, 15, 20, 25, 30, 40, 50, 60, 75, 100	50, 100	
120	0.1	0.02, 0.05, 0.1	5, 10, 15, 20, 25, 30, 40, 50, 60, 75, 100, 120	60, 120	
150	0.01, 0.02, 0.05, 0.1	0.05, 0.1	5, 10, 15, 20, 25, 30, 40, 50, 60, 75, 100, 120, 150	75, 150	
200	0.02, 0.01, 0.02, 0.05, 0.1, 0.2	0.05, 0.1, 0.2	5, 10, 15, 20, 25, 30, 40, 50, 60, 75, 100, 120, 150, 200	100, 200	
250	0.05, 0.1, 0.2	0.05, 0.1, 0.2	5, 10, 15, 20, 25, 30, 40, 50, 60, 75, 100, 120, 150, 200, 250	120, 250	
300	0.02, 0.05, 0.1, 0.2	0.05, 0.1, 0.2	5, 10, 15, 20, 25, 30, 40, 50, 60, 75, 100, 120, 150, 200, 250, 300	150, 300	
400	0.02, 0.05, 0.1, 0.2	0.1, 0.2	5, 10, 15, 20, 25, 30, 40, 50, 60, 75, 100, 120, 150, 200, 250, 300, 400	200, 400	
500	0.05, 0.1, 0.2, 0.5	0.1, 0.2, 0.5	5, 10, 15, 20, 25, 30, 40, 50, 60, 75, 100, 120, 150, 200, 250, 300, 400, 500	250, 500	
600	0.05, 0.1, 0.2, 0.5	0.1, 0.2, 0.5	5, 10, 15, 20, 25, 30, 40, 50, 60, 75, 100, 120, 150, 200, 250, 300, 400, 500, 600	300, 600	
750	0.05, 0.1, 0.2, 0.5	0.2, 0.5	5, 10, 15, 20, 25, 30, 40, 50, 60, 75, 100, 120, 150, 200, 250, 300, 400, 500, 600, 750	300, 750	
1000	0.05, 0.1, 0.2, 0.5, 1	0.2, 0.5, 1	5, 10, 15, 20, 25, 30, 40, 50, 60, 75, 100, 120, 150, 200, 250, 300, 400, 500, 600, 750, 1000	500, 1000	
1200	0.1, 0.2, 0.5, 1	0.2, 0.5, 1	5, 10, 15, 20, 25, 30, 40, 50, 60, 75, 100, 120, 150, 200, 250, 300, 400, 500, 600, 750, 1000, 1200	600, 1200	
1500	0.1, 0.2, 0.5, 1	0.5, 1	5, 10, 15, 20, 25, 30, 40, 50, 60, 75, 100, 120, 150, 200, 250, 300, 400, 500, 600, 750, 1000, 1200, 1500	750, 1500	
2000	0.1, 0.2, 0.5, 1, 2	0.5, 1, 2	5, 10, 15, 20, 25, 30, 40, 50, 60, 75, 100, 120, 150, 200, 250, 300, 400, 500, 600, 750, 1000, 1200, 1500, 2000	1000, 2000	
2500	0.2, 0.5, 1, 2	0.5 ,1, 2	5, 10, 15, 20, 25, 30, 40, 50, 60, 75, 100, 120, 150, 200, 250, 300, 400, 500, 600, 750, 1000, 1200, 1500, 2000, 2500	1200, 2500	
3000	0.2, 0.5, 1, 2	0.5 ,1 ,2	5, 10, 15, 20, 25, 30, 40, 50, 60, 75, 100, 120, 150, 200, 250, 300, 400, 500, 600, 750, 1000, 1200, 1500, 2000, 2500, 3000	1500, 3000	
5000	0.5, 1, 2, 5	1, 2, 5	5, 10, 15, 20, 25, 30, 40, 50, 60, 75, 100, 120, 150, 200, 250, 300, 400, 500, 600, 750, 1000, 1200, 1500, 2000, 2500, 3000, 5000	2500,5000	
6000	0.5, 1, 2, 5	1, 2, 5	5, 10, 15, 20, 25, 30, 40, 50, 60, 75, 100, 120, 150, 200, 250, 300, 400, 500, 600, 750, 1000, 1200, 1500, 2000, 2500, 3000, 5000, 6000	2500,5000	
7500	0.5, 1, 2, 5	2, 5	5, 10, 15, 20, 25, 30, 40, 50, 60, 75, 100, 120, 150, 200, 250, 300, 400, 500, 600, 750, 1000, 1200, 1500, 2000, 2500, 3000, 5000, 6000, 7500	3000,7500	
10000	0.5, 1, 2, 5, 10	2, 5, 10	5, 10, 15, 20, 25, 30, 40, 50, 60, 75, 100, 120, 150, 200, 250, 300, 400, 500, 600, 750, 1000, 1200, 1500, 2000, 2500, 3000, 5000, 6000, 7500, 10000	5000,10000	
12000	1, 2, 5, 10, 20	2, 5, 10	5, 10, 15, 20, 25, 30, 40, 50, 60, 75, 100, 120, 150, 200, 250, 300, 400, 500, 600, 750, 1000, 1200, 1500, 2000, 2500, 3000, 5000, 6000, 7500, 10000, 12000	6000,12000	
15000	1, 2, 5, 10	5, 10	5, 10, 15, 20, 25, 30, 40, 50, 60, 75, 100, 120, 150, 200, 250, 300, 400, 500, 600, 750, 1000, 1200, 1500, 2000, 2500, 3000, 5000, 6000, 7500, 10000, 12000, 15000	7500,15000	
20000	1, 2, 5, 10, 20	5, 10, 20	5, 10, 15, 20, 25, 30, 40, 50, 60, 75, 100, 120, 150, 200, 250, 300, 400, 500, 600, 750, 1000, 1200, 1500, 2000	10000,20000	

3.4.5 Graduation

Set the scale readability.

0.001, 0.002, 0.005, 0.01, 0.02, 0.05, 0.1, 0.2, 0.5, 1, 2, 5, 10, 20. **NOTE**: Not all settings are available for each capacity. Refer to the Setup Table 3.3 for available settings.

3.4.6 Power On Unit

Set the unit that will be active at power on.

oz, lb, g, kg, lb:oz or

Auto (last unit in use when power was turned off.)

3.4.7 Zero Range

Set the percentage of scale capacity that may be zeroed.

0% = zeroing disabled

2% = zero up to 2 percent of capacity

100% = zero up to full capacity

3.4.8 End Setup

Advance to the next menu.

3.5 Readout Menu

Enter this menu to customize display functionality.

Reset:	No, Yes	
Filter Level	Lo, Med,	Hi
Auto Zero Tracking	Off, 0.5d	, 1d, 3d
Backlight Auto Shut Off	Off, On, I	Auto
Auto Shut Off	Off	
End Readout	Exit READ	OUT menu

3.5.1 Reset

Set the Readout menu to factory default settings.

No = not reset

Yes = reset

If the Legal for Trade menu item is set to ON, the Stable Range, Averaging Level, Auto Zero Tracking and Auto Off settings are not reset.

Grßd
0.00 /
•
•
05

|--|

28-0	
<i>Û</i> -	8
<i>0</i> -	2
<i>[]</i> -	100

|--|

READOUT menu
r 8588

•	
	<i>no</i>
	465

3.5.2 Filter

Set the amount of signal filtering.

LO	0 = less stability, faster stabilization time (≤ 1 sec.)	
MEd	= normal stability, stabilization time (≤ 2 sec.)	

HI = greater stability, slower stabilization time (\leq 3 sec.)

3.5.3 Auto-Zero Tracking

Set the automatic zero tracking functionality.

	8
OFF	= disabled.
0.5 d	= the display will maintain zero until a drift of 0.5 divisions per second has been exceeded.
1 d	= the display will maintain zero until a drift of 1 division per second has been exceeded.
3 d	= the display will maintain zero until a drift of 3 divisions per second has been exceeded.

NOTE: When the LFT menu item is set to ON, the selections are limited to 0.5d and 3d. The setting is locked when the hardware lock switch is set to the ON position.

3.5.4 Backlight

Set the display backlight functionality.

. ,	
OFF	= always off.
ON	= always on.
AUtO	= turns on when a button is pressed or the displayed weight changes.
	turns off after 5 seconds of no activity.

3.5.5 Auto Off Timer

Set the automatic shut off functionality.

- OFF = disabled
- SEt 1 = powers off after 1 minute of no activity.
- SEt 2 = powers off after 2 minutes of no activity.
- SEt 5 = powers off after 5 minutes of no activity.



ÛFF	
0.5	ď
1	d
3	ď

L	'GHE
	OFF
	00
R	IUE D



3.5.6 End Readout

Advance to the next menu.

36001 3.6 Mode Menu Enter this menu to activate the desired application Reset: No, Yes modes. Off, On Count: End Mode Exit MODE menu

3.6.1 Reset Set the Mode menu to the factory defaults. rESEE No = not reset. Yes = reset. NOTE: If the Legal for trade menu item is set ON, the settings are not reset. YES 3.6.2 Parts Counting Mode Set the status. OFF = Disabled ON = Enabled

3.6.3 End Mode

Advance to the next menu.

COUNE OFF 00

00



3.7 Unit Menu

Enter this menu to activate the desired units. Default settings are bold.

	UN IF
Reset:	No, Yes
Kilograms:	Off, On
Pounds:	Off , On
Grams:	Off , On
Ounces:	Off , On
Pounds:Ounces	Off , On
End Unit	Exit UNIT menu

3.7.1 Reset

3.7.1 Reset	r E5EE
Set the Unit menu to the factory defaults.	
Settings:	
NO = not reset.	
YES = reset	985

If the Legal for Trade menu item is set ON, the settings are not reset.

3.7.2 Kilogram Unit

Set the status.

OFF	= Disabled
ON	= Enabled

3.7.3 Pound Unit

OFF	= Disabled
ON	= Enabled

3.7.4 Gram Unit

Set	the	status
001	1110	orarao.

OFF	= Disabled
ON	= Enabled

3.7.5 Ounce Unit

Set the status.

OFF = Disabled ON = Enabled

3.7.6 Pound Ounce Unit

Set the status.

OFF	= Disabled
ON	= Enabled

965
<i>∐∏ E</i> ₅g
000

OFF
00

UП IE '	b
0F <i>F</i>	
00	

UП IE "
OFF
00

ШΠ	! Ł ∞
01	c F
 l	0N

l	JП	lĿ	lb oz	
ÛFF				
	l	חכ		

3.7.7 End Unit

Advance to the next menu.

3.8 **Print Menu**

Pr int Enter this menu to define printing parameters. Default settings are bold. Reset No, Yes Baud Rate: 300, 600, 1200, 2400, 4800, 9600, 19200 Parity: 7 Even, 7 Odd, 7 None, 8 None 3.8.1 Reset rESEE Stop Bit 1 or 2 Set the Print menu to factory defaults. Handshake: Off, XON/XOFF = not reset. NO 00 Off, On Stable Only Auto Print Off, YES = reset. On Stable (-> Load, Load and Zero), *YES* Interval (-> 1...3600), Continuous Content Gross (->Off, On) NOTE: If the Legal for Trade menu item is set to ON, the following Net (->Off, On) settings are not reset: Stable, Auto Print Tare (->Off, On) Unit (->Off, On) End Print Exit PRINT menu

3.8.2 Baud

Set the Baud rate.


300	= 300 bps
600	= 600 bps
1200	=1200 bps
2400	= 2400 bps
4800	= 4800 bps
9600	= 9600 bps
19200	= 19200 bps

3.8.3 Parity

Set the data bits and parity.

7 EVEN = 7 data bits, even parity.7 Odd = 7 data bits, odd parity. 7 NONE = 7 data bits, no parity.8 NONE = 8 data bits, no parity.

6 <i>8</i> 87
300
600
1200
2400
4800
9600
19200
286 וצא
ח בטבח
7 079
3חסה ר
8 0002



End

3.8.4 Stop E Set the number of 1 2		SEOP 1 2
3.8.5 Hands Set the flow com NONE ON-OF	trol method. = no handshaking.	HRNJ NONE ON-OFF
3.8.6 Print Set the print crite OFF ON	Stable Data Only era. = values are printed immediately. = values are only printed when the stability criteria are met.	5ER6LE OFF 00
3.8.7 Auto F Set the automatia OFF ON.StAt INtEr CONt	c printing functionality. = disabled.	RPr int OFF ONSER5 INEEr CONE
	lected, set the Print Interval. 00 (seconds)	l 3600
3.8.8 Conte Select the addition GROSS NET TARE UNIT	onal content of the printout.	СОПЕПЕ GrOSS ЛЕЕ ЕЯгЕ ЦП IE

3.8.9 End Print

Advance to the next menu.

3.9 Menu	Lock Menu	<i>ב.רחבחט</i>
Enter this menu.	Default settings are bold.	Reset:No, YesLock Calibration MenuOff, OnLock Setup MenuOff, OnLock Readout MenuOff, OnLock Mode MenuOff, OnLock Unit MenuOff, OnLock Print MenuOff, OnEnd Lock MenuOff, On
3.9.1 Reset		r 8585
	sk menu to factory defaults.	
NO	= not reset.	00
YES	= reset.	985
NOTE: Settings for	or LFT controlled menu items are not reset.	
3.9.2 Lock C	Calibration	LEAL
Set the status.		
OFF ON	 Calibration menu is not locked. Calibration menu is locked and hidden. 	OFF
ÖN		00
3.9.3 Lock S Set the status.	Setup	L.SEEUP
OFF	= Setup menu is not locked.	OFF
ON	= Setup menu is locked and hidden.	
	Readout	L.r.E.R.d
Set the status. OFF	= Readout menu is not locked.	OFF
OFF	= Readout menu is locked and hidden.	
		00
3.9.5 Lock N	Node	1.000 <i>4E</i>
Set the status. OFF	= Mode menu is not locked.	ÛFF
ON	 Mode menu is locked and hidden. 	
		00
3.9.6 Lock U	Init	L.UN IE
Set the status. OFF	= Unit menu is not locked.	OFF
OFF	= Unit menu is locked and hidden.	
2		00

3.9.7 Lock Print

Set the status.

OFF	= Print menu is not locked.
ON	= Print menu is locked.

3.9.8 End Lock

Advance to the next menu.

3.10 Security Switch

A security switch is located on the Main PCB board. When the switch is set to the on position, user menu settings that were locked in the Menu Lock can not be changed.

Open the housing as explained in Section 2.3.1. Set the position of security switch to ON as shown in Figure 1-3.

4 **OPERATION**

4.1 Turning Indicator On/Off

To turn the Indicator on, press the and hold the **ON/ZERO** *Off* button for 2 seconds. The Indicator performs a display test, momentarily displays the software version, and then enters the active weighing mode.

To turn the Indicator off, press and hold the **ON/ZERO** Off button until OFF is displayed.

4.2 Zero Operation

Zero can be set under the following conditions:

- Automatically at Power On (initial zero).
- Semi-automatically (manually) by pressing the ON/ZERO Off button.
- Semi-automatically by sending the Zero command (Z or alternate zero command).


Press the **ON/ZERO** *Off* button to zero the weight display. The scale must be stable to accept zero operation.

4.3 Manual Tare

When weighing an item that must be held in a container, taring stores the container weight in memory. Place the empty container on the scale (example 0.5 kg) and press the **TARE** button. The display will show the net weight.

To clear the Tare value, empty the scale and press the **TARE** button. The display will show the gross weight.











	L.Pr	ιnΈ
	l]FF
		00
Γ		Ind

4.4 Changing Units of Measure

Press and hold the **PRINT** *Units* button until the desired measuring unit appears. Only measuring units enabled in the Unit Menu will be displayed (refer to Section 3.7).

4.5 Printing Data

Printing the displayed data to a printer or sending the data to a computer requires that the communication parameters in the Print Menu are set (refer to Section 3.8).

Press the **PRINT** *Units* button to send the displayed data to the communication port (the Auto-Print Mode in Section 3.8 function must be Off).

4.6 Application Modes

Only modes enabled in the mode menu will be displayed (refer to Section 3-6).

4.6.1 Weighing

Place the item to be weighed on the scale. The illustration indicates a sample of 1.5 kg, Gross weight.

Note: To return to the Weighing mode from the Parts Counting mode, press and hold the *Mode* button until WEIGH is displayed.

4.6.2 Parts Counting

Use this mode to count parts of uniform weight. The Indicator determines the quantity based on the average weight of a single part. All parts must be uniform in weight for accurate measurements.

To enter the Parts Counting mode, press and hold the *Mode* button until Count is displayed.

Average Piece Weight (APW)

When the *Mode* button is released, CLr.PW Pcs is displayed.

NOTE: If no APW has been previously stored, the CLr.PW display is skipped and the display shows PUt10Pcs.

Clearing a Stored APW

Press the Yes button to clear the stored APW.









Recalling a Stored APW

Press the No button to recall the existing APW.

Press the FUNCTION Mode button to temporarily display the APW value.





Establishing the Average Piece Weight (APW)

The display shows Put10 Pcs.

Establishing a New APW

Press the No button to increment the sample size. Choices are 5, 10, 20, 50, 100 and 200.

To establish the APW, place the specified quantity of samples on the scale and press the **Yes** button to capture the weight.

Begin Counting

Place the parts on the scale and read the count. If a container is used, be sure to tare the empty container first.

















5 SERIAL COMMUNICATION

The T31P and T31XW Indicators include an RS232 serial communication interface.

The setup of RS232 operating parameters are more fully explained in Section 3.8. The physical hardware connection is explained in in Section 2.2.

The interface enables display data to be sent to a computer or printer. A computer can be used to control some functions of the indicator using the commands listed in Table 5-1.

5.1 Interface Commands

Communicate to the indicator using the command characters listed in Table 5-1.

Command Character	Function				
	Immediate Drint of dianlayed weight (otable or unstable)				
IP	Immediate Print of displayed weight (stable or unstable).				
Р	Print stable displayed weight (according to stability setting).				
СР	Continuous Print.				
SP	Print when stable.				
хР	Interval Print x = Print Interval (1-3600 sec)				
Z	Same as pressing Zero button				
Т	Same as pressing Tare button				
хT	Download Tare value in grams (positive values only). Sending OT clears tare (if allowed)				
PU	Print current unit: g, kg, lb, oz, lb:oz				
хU	Set scale to unit x: 1=g, 2=kg, 3=lb, 4=oz, 5=lb:oz				
PV	Version: print name, software revision and LFT ON (if LFT is set ON).				
Esc R	Global reset to reset all menu settings to the original factory defaults				

TABLE 5-1. SERIAL INTERFACE COMMAND TABLE.

NOTES:

- Commands sent to the Indicator must be terminated with a carriage return (CR) or carriage return-line feed (CRLF).
- Data output by the Indicator is always terminated with a carriage return-line feed (CRLF).
- The xT (preset tare) command is not available when LFT is set to ON.

5.2 Output Format

The default serial output format is shown below.

Field:	Polarity	Space	Weight	Space	Unit	Stability	Legend	CR	LF
Length:	1	1	7	1	5	1	3	1	1

Definitions: Polarity, "-" sign if negative, blank if positive.

Weight, up to 6 numbers and 1 decimal, right justified, leading zero blanking. Units, up to 5 characters.

Stability, "?" character is printed if not stable, blank if stable.

Legend, up to 3 characters: G = gross weight, NET = net weight, T = tare

6. LEGAL FOR TRADE

When the indicator is used in trade or a legally controlled application it must be set up, verified and sealed in accordance with local weights and measures regulations. It is the responsibility of the purchaser to ensure that all pertinent legal requirements are met.

6.1 Settings

Before verification and sealing, perform the following steps:

- 1. Verify that the menu settings meet the local weights and measures regulations.
- 2. Perform a calibration.
- 3. Set Legal for Trade to ON in the Setup menu.
- 4. Exit the menu.
- 5. Disconnect power from the indicator and open the housing as explained in Section 2.3.1.
- 6. Set the position of the security switch to ON as shown in Section 1.2, Figure 1-3, Item 9.
- 7. Close the housing.
- 8. Reconnect power and turn the indicator on.

NOTE: For installations that employ the audit trail sealing method, steps 5 to 8 are not required. However, the security switch may be set to ON to safeguard against unintentional changes to configuration and calibration settings.

NOTE: When Legal for Trade is set to ON and the security switch is set to ON, the following menu settings cannot be changed: Span Calibration, Linearity Calibration, GEO, LFT, Calibration Unit, Capacity, Graduation, Power On Unit, Zero Range, Auto Zero Tracking, Expanded Mode, Count Mode, Kilogram Unit, Pound Unit, Gram Unit, Ounce Unit, Pound Ounce Unit, Stable Only. To enable editing of these menu settings, return the security switch to the off position and set LFT menu item to off.

6.2 Verification

The local weights and measures official or authorized service agent must perform the verification procedure. Please contact your local weights and measures office for further details.

6.3 Sealing

6.3.1 Physical Seals

For jurisdictions that use the physical sealing method, the local weights and measures official or authorized service agent must apply a security seal to prevent tampering with the settings. Refer to the illustrations below for sealing methods.




Figure 6-1. T31P Wire Seal




Figure 6-2. T31P Paper Seal





Figure 6-3. T31XW Wire Seal



6.3.2 Audit Trail Seal

For jurisdictions that use the audit trail sealing method, the local weights and measures official or authorized service agent must record the current configuration and calibration event counter values at the time of sealing. These values will be compared to values found during a future inspection.

NOTE: A change to an event counter value is equivalent to breaking a physical seal.

The audit trail uses two event counters to record changes to configuration and calibration settings.

- The configuration event counter (CFG) will index by 1 when exiting the menu if one or more of the following settings are changed Legal for Trade, Calibration Unit, Capacity, Graduation, Power On Unit, Zero Range, Auto Zero Tracking, Expanded Mode, Count Mode, Kilogram Unit, Pound Unit, Gram Unit, Ounce Unit, Pound Ounce Unit, Stable Only. Note that the counter only indexes once, even if several settings are changed. The configuration event counter values range from CFG000 to CFG999. When the value reaches CFG999, the count starts over at CFG000.
- The calibration event counter (CAL) will index by 1 when exiting the menu if a Span Calibration, Linearity Calibration or GEO setting change is made. Note that the counter only indexes once, even if several settings are changed. The calibration event counter values range from CAL000 to CAL999. When the value reaches CAL999, the count starts over at CAL000.

The event counters can be viewed by pressing and holding the MENU button. While the button is held, the display will show MENU followed by Audit.



Release the button when Audit is displayed to view the audit trail information.



3000 Series Indicators

The audit trail information is displayed in the format CFGxxx and CALxxx.





Then the indicator returns to normal operation.



7 MAINTENANCE



CAUTION: DISCONNECT THE UNIT FROM THE POWER SUPPLY BEFORE CLEANING.

7.1 Model T31P Cleaning

- The housing may be cleaned with a cloth dampened with a mild detergent if necessary.
- Do not use solvents, chemicals, alcohol, ammonia or abrasives to clean the housing or control panel.

7.2 Model T31XW Cleaning

- Use approved cleaning solutions for the stainless-steel Indicator housing and rinse with water. Dry thoroughly.
- Do not use solvents, chemicals, alcohol, ammonia or abrasives to clean the control panel.

7.3 Troubleshooting

SYMPTOM	PROBABLE CAUSE(s)	REMEDY
Unit will not turn on.	Power cord not plugged in or properly connected.	Check power cord connections. Make sure power cord is plugged in properly into the power outlet.
	Power outlet not supplying electricity.	Check power source.
	Battery power used up. (T31P Only)	Reconnect AC power to charge the battery.
	Other failure.	Service required.
Cannot zero the Scale, or will not zero when	Load on Scale exceeds allowable limits.	Remove load on Scale.
turned on.	Load on Scale is not stable.	Wait for load to become stable.
	Load Cell damage.	Service required.
Unable to calibrate.	Lock Calibration Menu set to On.	Set Lock Calibration Menu to Off. Refer to Section 3.9 Menu Lock.
	Lock switch is "on".	Set the Lock switch to Off.
	LFT menu set to On.	Set LFT menu to Off.
	Incorrect value for calibration mass.	Use correct calibration mass.
Cannot display weight in desired weighing unit.	Unit not set to On.	Enable unit in the Units Menu. Refer to Section 3.7 in the Unit Menu.
Cannot change menu settings.	Menu has been locked.	Set selected menu to Off in the Lock Menu. Lock Switch on the circuit board may need to be set to the Off position.
	Lock switch set on.	Set the Lock switch to off.
Battery indicator is flashing. (T31P Only)	Battery discharged.	Connect indicator to power and charge battery.
Battery fails to charge fully. (T31P Only)	Battery is defective.	Have the battery replaced by an authorized Ohaus service dealer.
Error 7.0	Unstable weight reading when defining reference weight.	Unstable Error, check platform location.

TABLE 7-1. TROUBLESHOOTING.

SYMPTOM	PROBABLE CAUSE(s)	REMEDY
Error 8.1	Weight reading exceeds Power On Zero limit.	Remove load from scale. Recalibrate scale.
Error 8.2	Weight reading below Power On Zero limit.	Add load to scale. Recalibrate scale.
Error 8.3	Weight reading exceeds Overload limit.	Reduce load on scale.
Error 8.4	Weight reading below Underload limit.	Add load to scale. Recalibrate scale.
Err 9.0	Internal fault	Service required.
Err 9.5	Calibration data not present.	Calibrate scale.
Err 53	EEPROM data incorrect.	Service required.
CAL E	Calibration Error. Calibration value outside allowable limits.	Repeat calibration using correct calibration weights.
LOW.rEF	The average piece weight of the parts is small (warning).	Use parts with average piece weight greater than or equal to 1 division.
REF.WT Err	The average piece weight of the parts is too small.	Use parts with a average piece weight greater than or equal to 0.1 division.

TABLE 7-1. TROUBLESHOOTING (Cont.).

7.4 Service Information

If the troubleshooting section does not resolve your problem, contact an authorized Ohaus Service Agent. For Service assistance in the United States, call toll-free 1-800-526-0659 between 8:00 AM and 5:00 PM Eastern Standard Time. An Ohaus Product Service Specialist will be available to assist you. Outside the USA, please visit our website www.ohaus.com to locate the Ohaus office nearest you.

8. TECHNICAL DATA

8.1 Specifications

Materials

T31XW Housing: stainless steel T31P Housing: ABS plastic Keypad: polyester Feet: Rubber Display Window: Polycarbonate

Ambient conditions

The technical data is valid under the following ambient conditions:

Ambient temperature:	-10°C to 40°C / 14°F to104°F	
Relative humidity:	Maximum relative humidity 80% for temperatures up to 31°C decreasing linearly to 50%	
	relative humidity at 40°C.	
Altitude:	up to 2000m	
Operability is assured at ambient temperatures between -10°C. and 40°C.		

TABLE 8-1. SPECIFICATIONS

Indicator	T31P	T31XW
Capacity Range	5 to 20000 kg or lb	
Maximum Displayed Resolution	1:20,000	
Type Approved Resolution	1:6,000	
Minimum Average Piece Weight (APW)	ld	
Weighing Units	kg, Ib, g,	oz, lb:oz
Functions	Weighing, Po	arts Counting
Display	1 in./2.5 cm digit height, 6-digit, 7-segment 1.5 in./3.8 cm high x 4.9 in./12.5 cm wide backlit LCD	
Backlight	White LED	
Keypad	4-button mechanical switches	4-button membrane switch
Ingress Protection		IP66
Load Cell Excitation Voltage	5V DC	
Load Cell Drive	Up to 4 x 350 ohm Load Cells	
Load Cell Input Sensitivity	Up to 3 mV/V	
Stabilization Time	Within 2 Seconds	
Auto-zero Tracking	Off, 0.5, 1 or 3 Divisions	
Zeroing Range	0%, 2% or 100% of Capacity	
Span Calibration	5 kg or 5 lb to 100% Capacity	
Interface	RS232	
Overall Dimensions (W x D x H) (in/mm)	8.2 x 2.8 x 6.5 / 210 x 71 x 168	8.3 x 2.8 x 5.8 / 212 x 71 x 149
Net Weight (Ib/kg)	3.6 / 1.6	6.6 / 2.9
Shipping Weight (Ib/kg)	5.7 / 2.6	8.8 / 4.0
Operating Temperature Range	-10°C to 40°C/14°F to 104°F	
Power	9 - 12VDC, 0.5A, AC Adapter with Internal rechargeable, Sealed Lead-Acid Battery (100-hour typical life) (T31P) 100-240 VAC / 50-60 Hz, Internal Power Supply (T31XW)	

8.2 Accessories

TABLE 8-2. ACCESSORIES.

DESCRIPTION	PART NUMBER
Column Mount Kit, 35 cm painted steel	80251743
Column Mount Kit, 70 cm painted steel	80251744
Column Mount Kit, 35 cm stainless steel	80251745
Column Mount Kit, 70 cm stainless steel	80251746
Wall Mount Kit, T31P	80251747
Wall Mount Kit, T31XW	80251748
Interface Cable/PC 25-pin, T31P	80500524
Interface Cable/PC 9-pin, T31P	80500525
Interface Cable/PC 9-pin, T31XW	80500552
Interface Cable/PC 25-pin, T31XW	80500553

8.3 Drawings and Dimensions




Figure 8-1. T31P Indicator Overall Dimensions.




Figure 8-2. T31XW Indicator Overall Dimensions with Mounting Bracket.

8.4 Compliance

Compliance

Compliance to the following standards is indicated by the corresponding mark on the product.

Marking	Standard	
<i>c c</i>	This product conforms to the EMC directive 2004/108/EC, the Low Voltage Directive 2006/95/EC and the	
	Non-automatic Weighing Instruments Directive 2009/23/EC. The complete Declaration of Conformity is	
	available online at www.ohaus.com.	
C	AS/NZS4251.1 Emission, AS/NZS4252.1 Immunity	
	UL60950-1: 2003	

EC Emissions Note

This device complies with EN55011/CISPR 11 Class B Group 1.

Important notice for verified weighing instruments

Weighing Instruments verified at the place of manufacture bear one of the preceding marks on the packing label and the green M' (metrology) sticker on the descriptive data plate. They may be put into service immediately.



Weighing Instruments to be verified in two stages have no green 'M' (metrology) on the descriptive data plate and bear one of the preceding identification marks on the packing label. The second stage of the initial verification must be carried out by an authorized and certified service organization established within the European Community or by the National Notified Body.

The first stage of the initial verification has been carried out at the manufacturer's work. It comprises all tests according to the adopted European Standard EN 45501:1992, paragraph 8.2.2.

If national regulations limit the validity period of the verification, the user of the weighing instrument must strictly observe the re-verification period and inform the respective Weights and Measures authority.



Disposal

In conformance with the European Directive 2002/96/EC on Waste Electrical and Electronic Equipment (WEEE) this device may not be disposed of in domestic waste. This also applies to countries outside the EU, per their specific requirements.

The Batteries Directive 2006/66/EC introduces new requirements from September 2008 on removability of batteries from waste equipment in EU Member States. To comply with this Directive, this device has been designed for safe removal of the batteries at end-of-life by a waste treatment facility.

Please dispose of this product in accordance with local regulations at the collecting point specified for electrical and electronic equipment.

If you have any questions, please contact the responsible authority or the distributor from which you purchased this device.

Should this device be passed on to other parties (for private or professional use), the content of this regulation must also be related.

For disposal instructions in Europe, refer to www.ohaus.com, choose your country then search for WEEE.

Thank you for your contribution to environmental protection.

FCC Note

This equipment has been tested and found to comply with the limits for a Class A digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference in which case the user will be required to correct the interference at his own expense.

Industry Canada Note

This Class A digital apparatus complies with Canadian ICES-003.

ISO 9001 Registration

In 1994, Ohaus Corporation, USA, was awarded a certificate of registration to ISO 9001 by Bureau Veritus Quality International (BVQI), confirming that the Ohaus quality management system is compliant with the ISO 9001 standard's requirements. On May 21, 2009, Ohaus Corporation, USA, was re-registered to the ISO 9001:2008 standard.

LIMITED WARRANTY

Ohaus products are warranted against defects in materials and workmanship from the date of delivery through the duration of the warranty period. During the warranty period Ohaus will repair, or, at its option, replace any component(s) that proves to be defective at no charge, provided that the product is returned, freight prepaid, to Ohaus.

This warranty does not apply if the product has been damaged by accident or misuse, exposed to radioactive or corrosive materials, has foreign material penetrating to the inside of the product, or as a result of service or modification by other than Ohaus. In lieu of a properly returned warranty registration card, the warranty period shall begin on the date of shipment to the authorized dealer. No other express or implied warranty is given by Ohaus Corporation. Ohaus Corporation shall not be liable for any consequential damages.

As warranty legislation differs from state to state and country to country, please contact Ohaus or your local Ohaus dealer for further details.



Ohaus Corporation 7 Campus Drive Suite 310 Parsippany, NJ 07054, USA Tel: (973) 377-9000 Fax: (973) 944-7177 www.ohaus.com



P/N 80251064 E © 2012 Ohaus Corporation, all rights reserved.

Printed in China



Defender™ T Series Bases Instruction Manual



Defender Series Base

INTRODUCTION

This manual covers installation and maintenance instructions for the Ohaus Defender[™] T Series Bases. Please read this manual completely before installation and operation.

SAFETY PRECAUTIONS

For safe and dependable operation of this product, please comply with the following precautions:



- Use caution when lifting or moving heavy objects.
- Operate the base only under ambient conditions specified in this manual.
- Engage wheel brakes before operating the base.
- Ensure that the load cell cable does not pose an obstruction or tripping hazard.
- Do not operate the base in hazardous environments or unstable locations.
- Do not drop loads on the base.
- Do not lift the base by the top frame. Always lift from the bottom frame when moving the base.
- Service should only be performed by qualified personnel.



FAILURE TO FOLLOW THESE PRECAUTIONS COULD RESULT IN DAMAGE TO THE EQUIPMENT OR BODILY HARM.

INSTALLATION

Unpacking

Unpack and inspect the product to make sure that all components have been included. The package includes the following:

- Base
- Weighing Pan
- Instruction Manual

When purchased as a complete Defender Series Scale, the package will also include:

- Indicator
- Column Assembly

Note: Remove the red plastic shipping tabs from the corners of the scale base.

Assembly

Weighing Pan

Place the weighing pan securely over the rubber load pads on the top frame of the base.

Wiring Connections

Install the load cell cable to an indicator using the wiring codes in Table 1.

Note: When purchased as a Defender Series Scale, the base is already pre-wired to the indicator.

TABLE T. LUAD GELL CONNECTION	
FUNCTION	WIRE COLOR
+ Excitation	Green
 Excitation 	Black
+ Signal	White
– Signal	Red
Ground	Yellow

TABLE 1. LOAD CELL CONNECTION

Selecting the Location

To ensure accuracy, proper performance and safety, locate and operate the base on a stable, level surface. Avoid locations with rapid temperature changes or excessive dust, air currents, vibrations, electromagnetic fields or heat.

Locate the base on a level surface, making sure that the bubble in the level indicator is centered. When the base is level, engage the wheel brakes of the two rear wheels.

Note: Ensure that the base is level each time its location is changed.



MAINTENANCE

Cleaning

The base components should be kept clean and free of excessive material build up. A damp cloth with water and a mild detergent may be used to wipe clean the external surfaces – do not use acids, alkalis, strong solvents or abrasive materials.

Troubleshooting

Operational difficulties that may be encountered can often be traced to simple causes such as:

- Loose or incorrect wiring connections
- Obstructions to the base frame
- Unstable environments
- Incorrect calibration or setup of the indicator

If the troubleshooting section does not resolve or describe your problem, contact your authorized Ohaus service agent. For service assistance or technical support in the United States call toll-free 1-800-526-0659 between 8:00 AM and 5:00 PM EST. An Ohaus product service specialist will be available to provide assistance. Outside the USA, please visit our web site, www.ohaus.com to locate the Ohaus office nearest you.

TECHNICAL DATA

Technical Specifications

The technical data is valid under the following ambient conditions: Ambient temperature: -10° C to $+40^{\circ}$ C Relative humidity: 10% to 95% relative humidity, non-condensing Height above sea level: Up to 2,000m

	TABLE 2. SPECIFICATIONS	
MODEL	D250TX	D500TX
Capacity	250 kg / 500 lb	500 kg / 1000 lb
Approved Resolution	2500d	
Non-approved Resolution	5000d	
Safe Overload Capacity	150% of capacity	
Pan Dimensions	508 x 710 mm / 20 x 28 in.	
Base Dimensions	520 x 720 x 198 mm / 20.5 x 28.3 x 7.8 in.	
Base Construction	304 stainless steel platform with painted steel frame and polyurethane wheels	
Repeatability (std. deviation)	ld	
Linearity	±ld	
Load Cell Capacity	500 kg	1000 kg
Load Cell Cable	2 m L x 4-wire	
Load Cell Type	350 Ohm, aluminum, single point	
Load Cell Excitation	5-15V DC/AC	
Load Cell Rated Output	2mV/V	
Net Weight	40 kg / 88 lb	
Shipping Weight	46 kg / 101 lb	

Drawing




Figure 1. Defender Base Dimension Drawing.

Compliance

Compliance to the following standards is indicated by the corresponding marking on the product.

Marking	Standard
CE	This product conforms to the EMC directive 2004/108/EC and the Low Voltage Directive 2006/95/EC. The complete declaration of Conformity is available from Ohaus Corporation.
CC 09-056	NIST Handbook 44
AM-5723	Weights and Measures Regulations of Canada



Disposal

In conformance with the European Directive 2002/96 EC on Waste Electrical and Electronic Equipment (WEEE) this device may not be disposed of in domestic waste. This also applies to countries outside the EU, per their specific requirements.

Please dispose of this product in accordance with local regulations at the collecting point specified for electrical and electronic equipment.

If you have any questions, please contact the responsible authority or the distributor from which you purchased this device.

Should this device be passed on to other parties (for private or professional use), the content of this regulation must also be related.

For disposal instructions in Europe, refer to www.ohaus.com/weee.

Thank you for your contribution to environmental protection.

ISO 9001 Registration

In 1994, Ohaus Corporation, USA, was awarded a certificate of registration to ISO 9001 by Bureau Veritus Quality International (BVQI), confirming that the Ohaus quality management system is compliant with the ISO 9001 standard's requirements. On May 15, 2003, Ohaus Corporation, USA, was re-registered to the ISO 9001:2000 standard.

LIMITED WARRANTY

Ohaus products are warranted against defects in materials and workmanship from the date of delivery through the duration of the warranty period. During the warranty period Ohaus will repair, or, at its option, replace any component(s) that proves to be defective at no charge, provided that the product is returned, freight prepaid, to Ohaus.

This warranty does not apply if the product has been damaged by accident or misuse, exposed to radioactive or corrosive materials, has foreign material penetrating to the inside of the product, or as a result of service or modification by other than Ohaus. In lieu of a properly returned warranty registration card, the warranty period shall begin on the date of shipment to the authorized dealer. No other express or implied warranty is given by Ohaus Corporation. Ohaus Corporation shall not be liable for any consequential damages.

As warranty legislation differs from state to state and country to country, please contact Ohaus or your local Ohaus dealer for further details.



Ohaus Corporation 19A Chapin Road P.O. Box 2033 Pine Brook, NJ 07058-2033, USA Tel: (973) 377-9000 Fax: (973) 944-7177

With offices worldwide. www.ohaus.com



P/N 80253148 © 2009 Ohaus Corporation, all rights reserved.

Printed in China