Operating Instructions

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
Mode d'emploi

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Istruzioni d'uso

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FiveEasy™ FiveGo™
Conductivity meter FE30 / FG3

METTLER TOLEDO

English

Deutsch

Français

Español

Italiano

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Contents

1.	Introduction	3
2.	Safety measures	3
3.	Installation	4
3.1	Installing the electrode arm (optional FE30)	4
3.2	Installing the batteries (FG3)	4
3.3	FiveGo™ clip (FG3)	5
3.4	Installing the IP54 cover (FG3)	5
4.	Operating the FiveEasy™ and FiveGo™ conductivity meter	6
4.1	Display and key controls	6
4.2	Calibration	8
4.2.1	Selecting a standard	8
4.2.2	Performing a calibration	8
4.3	Sample measurement	8
4.4	TDS / salinity measurement	9
4.5	Setup	9
4.5.1	Set temperature correction coefficient	9
4.5.2	Set reference temperature	9
4.5.3	Set temperature unit	9
4.5.4	Set TDS factor	9
4.6	Using the memory (FG3 only)	10
4.6.1	Storing a reading	10
4.6.2	Recalling from memory	10
4.6.3	Clearing the memory	10
4.7	Self-diagnosis	10
4.8	Error message	11
5.	Maintenance	11
5.1	Meter maintenance	11
5.2	Disposal	11
6.	Accessories	12
7.	Specifications	13
8.	Appendix	15
8.1	Conductivity standards	15
8.2	Examples of temperature coefficients (α -value)	15
8.3	Practical salinity scale (UNESCO 1978)	16
8.4	Conductivity to TDS conversion factors	16
8.5	Error Limits	17

1. Introduction

Thank you for purchasing this high quality METTLER TOLEDO meter. Ease of use, reliable results and ergonomic design – these are our goals for products.

The Five meters have an excellent price/performance ratio and many useful accessories (e.g. electrode clip and wrist strap for the portable meter) are already included in the standard delivery of the instruments. Other accessories such as electrodes, solutions, an electrode arm for the benchmeter, or a carry bag for the portable meter are of course also available from METTLER TOLEDO.

You will find much helpful functionality in these meters: Automatic endpoint format which freezes the measurement automatically when the value is stable, predefined conductivity standards and IP54 design (FG3) are only few of many exiting features.

2. Safety measures

Measures for your protection

- Never work in an environment subject to explosion hazards! The housing of the instrument is not gas tight (explosion hazard due to spark formation, corrosion caused by the ingress of gases).

- When using chemicals and solvents, comply with the instructions of the producer and the general lab safety rules!

Measures for operational safety



- Do not unscrew the two halves of the housing!
- Dry off any liquid spills immediately! The instrument is not watertight (FG3 is IP54)!
- For FG3 instruments use batteries of the specified type only. Otherwise, proper operation cannot be guaranteed.
- Exclude the following environmental influences:
 - Powerful vibrations,
 - Direct sunlight,
 - Atmospheric humidity greater than 80%,
 - Corrosive gases,
 - Temperatures below 5 °C and above 40 °C,
 - Powerful electric or magnetic fields!

3. Installation




Carefully unpack the meter. Keep the calibration certificate in a safe place.

Insert the right adapter clip into the power adapter slot:



3.1 Installing the electrode arm (optional FE30)

The electrode arm can be fastened to FiveEasy™ on the left or the right side.


1. Remove the electrode stand rubber cover on that side of the bench meter where you prefer the stand.
Keep the cover in a safe place;
2. Insert the electrode stand firmly into the opening as shown;
3. Fix the opening of the arm upon the stand from above, adjust it to a moderate height.

3.2 Installing the batteries (FG3)

1. Pull the fixing point of the battery cover in the direction of the arrow, hold the lid with two fingers and remove it;
2. Insert the batteries in the battery compartment as shown;
3. Replace the battery cover.


3.3 FiveGo™ clip (FG3)

The FiveGo™ clip (Part no.51302960) is an electrode holder that can be placed next to the display on either side of the housing. To mount the clip, remove the cover from the clip's fixing point. Attach the clip by pressing it into the recess. Slide the shaft of the sensor into the clip from above.


3.4 Installing the IP54 cover (FG3)

Use the small plastic tool as shown:

4. Operating the FiveEasy™ and FiveGo™ conductivity meter

4.1 Display and key controls


- 1 Calibration icon - calibration in progress
- 2 Measurement icon - measurement or calibration is running
- 3 Setup icon - instrument is in the setup mode
- 4 Battery status icon (FG3 only) shows the condition of the batteries - fully charged, half-charged or fully discharged (To replace batteries, see section 3.2)
- 5 TDS factor in the setup mode
- 6 Conductivity / TDS / Salinity reading / Cell constant in calibration process
- 7 Endpoint stability / Auto endpoint
- 8 Temperature during measurement or reference temperature in setup mode
- 9 Standard  / Memory number **MR** (FG3 only) / Error index **Err**

Keys for FE30

	Press & release 	Press & hold for 3 seconds 
	- Start or endpoint measurement - Confirm setting, store entered value	- Turn auto endpoint on / off $\sqrt{A}$, $\sqrt{\quad}$
	- Start calibration	- Review the latest calibration data
	- Meter on - Back to measurement screen	- Meter off
	- Enter setup mode - Increase value during setting	
	- Switch between conductivity, TDS and salinity measuring modes - Decrease value during setting	
	- Start self-diagnosis	

Keys for FG3

	Press & release 	Press & hold for 3 seconds 
	- Start or endpoint measurement - Confirm setting, store entered value	- Turn auto endpoint on / off $\sqrt{A}$, $\sqrt{\quad}$
	- Start calibration	- Review the latest calibration data
	- Meter on - Back to measurement screen	- Meter off
	- Store current reading to memory - Increase value during setting - Scroll up through the memory	- Recall stored data
	- Switch between conductivity, TDS and salinity measuring modes - Decrease value during setting - Scroll down through the memory	- Enter setup mode
	- Start self-diagnosis	

4.2 Calibration

4.2.1 Selecting a standard

When using the Five conductivity meter, you have to select a standard for calibration.

Press the **Setup** key (FE30) or press and hold the **Mode/Setup** key (FG3) until the setup icon appears on the display and the current standard blinks. Use **▲** or **▼** to select another standard and press **Read** to confirm your selection. Press **Exit** to leave the setup mode.

The three predefined standards are:

- 84 $\mu\text{S/cm}$
- 1413 $\mu\text{S/cm}$
- 12.88 mS/cm

Tables for automatic temperature compensation are programmed in the meter for each standard (see also appendix).

4.2.2 Performing a calibration

Place the conductivity sensor in the defined calibration standard and press **Cal**.

The calibration icon and the measurement icon appear on the display. The meter endpoints according to the preselected endpoint mode after the signal has stabilized or after pressing **Read**. The standard value at the measurement temperature is displayed and stored; the measurement icon disappears from the display.

The cell constant is shown on the display for 3 seconds. To reject the calibration, press **Exit**, before the meter returns automatically after 3 seconds to the measurement screen.

Note:

To ensure the most accurate conductivity readings, you should verify your cell constant with a standard solution once a day and recalibrate if necessary. Use always fresh standards.

4.3 Sample measurement

Place the conductivity sensor in the sample and press **Read** to start the measurement: the measurement icon appears on the display and the decimal point blinks. The display shows the conductivity of the sample. The automatic endpoint **A** is the default setting of the meter. When the signal has stabilized, the display freezes automatically, and **√A** appears.

By pressing and holding **Read**, you can switch between the Auto and Manual endpoint mode. To manually endpoint a measurement, press **Read**: the display freezes and **√** appears.

Stability criterion for conductivity measurement – The sensor input signal of the meter may not change by more than 0.4% from the measured average conductivity of the probe in 6 seconds.

4.4 TDS / salinity measurement

To perform a TDS / salinity measurement, follow the same procedure as for a conductivity measurement. Press **Mode** to switch between conductivity, TDS and salinity measurement mode.

4.5 Setup

4.5.1 Set temperature correction coefficient

Press the **Setup** key (FE30) or press and hold the **Mode/Setup** key (FG3) until the setup icon appears on the display and the current standard blinks. Press **Read** to ignore standard setting; the meter changes to the α -coefficient setting automatically.

Use **▲** or **▼** to increase or decrease the temperature correction value. Press **Read** to confirm your setting. Continue with reference temperature setting or press **Exit** to go back to the measurement screen.

4.5.2 Set reference temperature

After confirming the temperature correction coefficient, the reference temperature appears on the display. Use **▲** or **▼** to toggle between 25 °C and 20 °C. Press **Read** to confirm your selection. Continue with temperature unit setting or press **Exit** to go back to the measurement screen.

4.5.3 Set temperature unit

After selecting the reference temperature, the current temperature unit blinks. Use **▲** or **▼** to toggle between °C and °F. Press **Read** to confirm your selection. Continue with TDS factor setting or press **Exit** to go back to the measurement screen.

4.5.4 Set TDS factor

After selecting the temperature unit, the current TDS factor blinks. Use **▲** or **▼** to increase or decrease the value. Press **Read** to confirm your setting. The meter will exit automatically to the measurement screen.

4.6 Using the memory (FG3 only)

4.6.1 Storing a reading

The FiveGo™ FG3 conductivity meter can store up to 30 measured results. Press **STO** when the measurement has endpointed. **M01** indicates that one result has been stored. If you press **STO** when **M30** is displayed, **FUL** indicates the memory is full. To store further data you will have to clear the memory (see below.)

4.6.2 Recalling from memory

Press and hold **RCL** to recall the stored values from the memory when the current measurement has endpointed.

Press **▲** or **▼** to scroll through the stored results. **R01** to **R30** indicates which result is being displayed. Press **Read** to exit.

4.6.3 Clearing the memory

Continue pressing **▲** or **▼** to scroll through the stored results until **MRCL** appears. Then press **Read**, **CLr** blinks. Press **Read** again to confirm the deletion or press **Exit** to return to the measurement mode without deleting the data.

4.7 Self-diagnosis

Press and hold **Read** and **Cal** simultaneously until the meter displays the full screen. Each icon blinks one after the other. This way you may check whether all icons are correctly shown. The next step is to check that the keys are functioning correctly. This requires user interaction.

When **b** blinks, five icons are displayed. Press the five keys in any order. Each time you press a key an icon disappears from the screen, continue to press the other keys until all the icons have disappeared.

When the self-diagnosis has been completed successfully, **PAS** appears. If self-diagnosis fails, error message **Err 1** appears (see section 4.8).

Note:

You have to finish pressing all five keys within two minutes, otherwise **Err 1** appears and you will have to repeat the procedure.

4.8 Error message

Error 0	Memory access error	Reset to factory settings
Error 1	Self-diagnosis failed	Repeat the self-diagnosis procedure and make sure that you finish pressing all five keys within two minutes.
Error 2	Measured values out of range	Check if the electrode is properly connected and placed in a sample solution.
Error 3	Measured standard temperature out of range (5 ... 35 °C)	Keep the standard temperature within the range for calibration.
Error 4	Measuring temperature out of range (0 ... 100 °C)	Check if the electrode is properly connected and keep the sample temperature within the range.
Error 9	The current data set has already been stored once	A measurement can only be stored once. Perform a new measurement to store a new data set.

5. Maintenance

5.1 Meter maintenance

Never unscrew the two halves of the housing!

The Five series instruments do not require any maintenance other than an occasional wipe with a damp cloth and the replacement of depleted batteries (FG3).

The housing is made of acrylonitrile butadiene styrene/polycarbonate (ABS/PC). This material is attacked by some organic solvents, such as toluene, xylene and methyl ethyl ketone (MEK). Any spillage should be immediately wiped off.

5.2 Disposal

In conformance with the European Directive 2002/96/EC on Waste Electrical and Electronic Equipment (WEEE) this device may not be disposed of in domestic waste. This also applies to countries outside the EU, per their specific requirements.

Please dispose of this product in accordance with local regulations at the collecting point specified for electrical and electronic equipment.

If you have any questions, please contact the responsible authority or the distributor from which you purchased this device.

Should this device be passed on to other parties (for private or professional use), the content of this regulation must also be related.

Thank you for your contribution to environment protection.

6. Accessories

Item	Order No.
FiveEasy™ conductivity meter	51302912
FiveEasy™ Power supply	51302950
FiveEasy™ electrode arm	51302951
Rubber covers (for electrode arm hole)	51302952
FiveGo™ conductivity meter	51302903
FiveGo™ electrode clip	51302960
Rubber clip covers (for electrode clip slot)	51302961
FiveGo™ Battery cover	51302962
FiveGo™ Wrist strap	51302964
FiveGo™ carry bag	51302966
IP54 electrode caps	51302965
Sealing kit (O-rings for battery case)	51302967
Quick guide FE30	51710493
Quick guide FG3	51710495
84 µS/cm standard solution, 500mL	51302153
1413 µS/cm standard solution sachets, 30 x 20 mL	51302049
1413 µS/cm standard solution, 250 mL	51300138
12.88 mS/cm standard solution sachets, 30 x 20 mL	51302050
12.88 mS/cm standard solution, 250 mL	51300139
LE703, conductivity sensor	51340335

7. Specifications

	FiveEasy™ conductivity meter FE30	FiveGo™ conductivity meter FG3
Measurement range		
Conductivity	0.0 µS/cm...199.9 mS/cm	0.0 µS/cm...199.9 mS/cm
TDS	0.1 mg/L ... 199.9 g/L	0.1 mg/L ... 199.9 g/L
Salinity	0.00 ... 19.99 psu	0.00 ... 19.99 psu
Temperature	0° C ... 100 °C	0° C ... 100 °C
Resolution		
Conductivity	Auto range	Auto range
	0.10 µS/cm ... 19.99 µS/cm	0.10 µS/cm ... 19.99 µS/cm
	20.0 µS/cm ... 199.9 µS/cm	20.0 µS/cm ... 199.9 µS/cm
	200 µS/cm ... 1999 µS/cm	200 µS/cm ... 1999 µS/cm
	2.00 mS/cm ... 19.99 mS/cm	2.00 mS/cm ... 19.99 mS/cm
	20.0 mS/cm ... 199.9 mS/cm	20.0 mS/cm ... 199.9 mS/cm
TDS	Auto range, same as conductivity	Auto range, same as conductivity
Salinity	0.00 ... 19.99 psu	0.00 ... 19.99 psu
Temperature	0.1 °C	0.1 °C
Limits of error		
Conductivity	± 0.5 % of measured value	± 0.5 % of measured value
TDS	± 0.5 % of measured value	± 0.5 % of measured value
Salinity	± 0.5 % of measured value	± 0.5 % of measured value
Temperature	± 0.3 °C	± 0.3 °C
Power requirements	9 V DC, 0.4 W ± 10% For use with CSA Certified (or equivalent approved) power source, which must have a limited circuit output.	Rating: 6VDC, 5mA Powered with batteries: 4xAAA/LR3 1.5V or NiMH 1.2 rechargeable
Size/weight	200 x 175 x 52 mm / 0.6 kg	169 x 82 x 36 mm / 0.18 kg (without batteries)
Display	Liquid crystal	Liquid crystal
Signal input	Mini-Din	Mini-Din

Specifications

IP rating	For indoor use only.	IP54 for system
Battery life	-	>250 working hours
Ambient conditions	Temperature: 5...40 °C	Temperature: 5...40 °C
	Relative humidity: 5 %...80 % (non-condensing)	Relative humidity: 5 %...80 % (non-condensing)
	Height up to 2000 m	Height up to 2000 m
	Installation category: II	Installation category: II
	Pollution degree: 2	Pollution degree: 2
Materials	Housing: ABS/PC reinforced	Housing: ABS/PC reinforced
	Window: Membrane	Window: polymethylmethacrylate (PMMA)
	Keypad: Membrane	Keypad: Membrane

8. Appendix

8.1 Conductivity standards

T(°C)	84 µS/cm	1413 µS/cm	12.88 mS/cm
5	53 µS/cm	896 µS/cm	8.22 mS/cm
10	60 µS/cm	1020 µS/cm	9.33 mS/cm
15	68 µS/cm	1147 µS/cm	10.48 mS/cm
20	76 µS/cm	1278 µS/cm	11.67 mS/cm
25	84 µS/cm	1413 µS/cm	12.88 mS/cm
30	92 µS/cm	1552 µS/cm	14.12 mS/cm
35	101 µS/cm	1667 µS/cm	15.39 mS/cm

8.2 Examples of temperature coefficients (α -value)

Substance at 25 °C	Concentration [%]	Temperature coefficient α [%/°C]
HCl	10	1.56
KCl	10	1.88
CH ₃ COOH	10	1.69
NaCl	10	2.14
H ₂ SO ₄	10	1.28
HF	1.5	7.20

α -coefficients of conductivity standards for a calculation to reference temperature of 25°C

Standard	Measurement temp.: 15°C	Measurement temp.: 20°C	Measurement temp.: 30°C	Measurement temp.: 35°C
84 µS/cm	1.95	1.95	1.95	2.01
1413 µS/cm	1.94	1.94	1.94	1.99
12.88 mS/cm	1.90	1.89	1.91	1.95

8.3 Practical salinity scale (UNESCO 1978)

In the Five conductivity meters the salinity is calculated according to the official definition of UNESCO 1978. Therefore the salinity S_{psu} of a sample in psu (practical salinity unit) at standard atmospheric pressure is calculated as follows:

$$S = \sum_{j=0}^5 a_j R_T^{j/2} - \frac{(T-15)}{1+k(T-15)} \sum_{j=0}^5 b_j R_T^{j/2}$$

$a_0 = 0.0080$	$b_0 = 0.0005$	$k = 0.00162$
$a_1 = -0.1692$	$b_1 = -0.0056$	
$a_2 = 25.3851$	$b_2 = -0.0066$	
$a_3 = 14.0941$	$b_3 = -0.0375$	
$a_4 = -7.0261$	$b_4 = 0.0636$	
$a_5 = 2.7081$	$b_5 = -0.0144$	

$$R_T = \frac{C_{\text{sample}}(T)}{C_{\text{KCl}}(T)}$$

(32.4356 g KCl per 1000 g of solution)

8.4 Conductivity to TDS conversion factors

Conductivity At 25 °C	TDS KCl		TDS NaCl	
	ppm value	Factor	ppm value	Factor
84 µS	40.38	0.5048	38.04	0.4755
447 µS	225.6	0.5047	215.5	0.4822
1413 µS	744.7	0.527	702.1	0.4969
1500 µS	757.1	0.5047	737.1	0.4914
8974 µS	5101	0.5685	4487	0.5000
12.880 µS	7447	0.5782	7230	0.5613
15.000 µS	8759	0.5839	8532	0.5688
80 mS	52.168	0.6521	48.384	0.6048

8.5 Error Limits

Message	Description	Range not accepted
Err 2	Measured values out of range	C: > 199.9 mS/cm TDS: < 0.1 mg/L or > 199.9 g/L SAL: > 19.99 psu
Err 3	Calibration standard temperature out of range	T: < 5 °C or > 35 °C
Err 4	Temperature out of range	T: < 0 °C or > 100 °C

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Quality certificate. Development, production and testing according to ISO9001.
Environmental management system according to ISO14001.

SERVICE

Worldwide service. Our extensive service network is among the best in the world and ensures maximum availability and service life of your product.

European conformity. The CE conformity mark provides you with the assurance that our products comply with the most recent EU directives.

INTERNET

On the Internet. You will quickly find lots of essential information about our products, our services, and our company at
<http://www.mt.com>

P51710489

Subject to technical changes and to the availability
of the accessories supplied with the instruments.

© Mettler-Toledo AG 2007 ME-51710489 Printed in China 0702/2.52

Mettler-Toledo AG, Analytical, Sonnenbergstrasse 74, CH-8603 Schwerzenbach,
Tel. +41 44 806 77 11, Fax +41 44 806 73 50, Internet:<http://www.mt.com>