**General instruction for pH electrodes**  
**Instrucciones generales para electrodos de pH**  
**Instructions générales pour**  
**les électrodes de pH**

Ohaus Corporation  
7 Campus Drive  
Suite 310  
Parsippany, NJ 07054 USA  
Tel: (973) 377-9000  
Fax: (973) 944-7177

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**Specifications:****Connection:** BNC**reference system:**

double salt bridge, Ag/AgCl--KCl

**Salt bridge solution:**

saturated Pure KCl solution

**Cable length:** 1.0m**Shaft length:** 120mm**Shaft diameter:** 12mm**pH range:** 0-14pH**Zero potential:**  $7.00 \pm 0.25$ pH**Slope:**  $\geq 97\%$ **Impedance:**  $\leq 250$ M $\Omega$ **Shaft material:** Glass**Reference Electrolyte:**

3M KCl Solution

**Temperature range:** 0-100°C**Reference Diaphragm:** Sleeve

**OHAUS** offer different pH electrode for different application. Make sure you have chosen the best fit pH electrode for your samples; and do calibration first for the electrode.

- 2-in-1 pH electrodes need to work with STTEMP30 for 0.01 pH accuracy.
- ST260 is fit for the sample which may react with Ag<sup>+</sup> ions such as tri-buffer solution.

**Before** using electrode in calibration or measurement, we should **note**:

1. When unpacking, check the electrode for any mechanical damage to the membrane, diaphragm or body.
2. Screw off watering cap first, then the lid, rinse electrode with pure water. Wipe electrode with a tissue. Do not rub the glass membrane, or the response time might increase.
3. Eliminate any air bubbles inside the membrane and the salt bridge solution. Make the electrode inside pressure equilibrium with the atmosphere by opening the electrode refilling aperture when measurement.

### **FAQs:**

#### **1. What's the meaning of electrode condition icon?**

After calibration, the pH meter will judge the electrode slope and offset, then display the judgment in the screen:

😊 slope > 95%

Electrode condition is very good

😐 slope : 90%-95%

Electrode condition is acceptable

😞 slope : 85%-90%

Electrode condition is not good, need re-calibration or replace

**2. How do you properly store the pH electrode after every use?**

A: Rinse the electrode with pure water, put the electrode into the watering cap. If the solution in the cap is polluted or dry, clean it and fill it with 3M KCl solution.

**3. What's the life of the pH electrodes?**

A: Measured samples, storage conditions and maintenance methods all could influence the life of pH electrodes; normally a pH electrode expires after 12 months. Some sample, such as strong acid or alkali, corrosive liquid will shorten the pH electrode life. pH electrode is consumable, is generally provided warranty.

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Ohaus Corporation, all rights reserved/ todos