


ORP half-cell



Order number	59904408
---------------------	----------

former order number	108503079
---------------------	-----------

Specifications

pH range	-
----------	---

zero point (pH)	-
-----------------	---

temperature range	-30 - 80 °C
-------------------	-------------

temperature probe	-
-------------------	---

shaft material	polypropylene
----------------	---------------

membrane glass	-
----------------	---

membrane resistance	-
---------------------	---

reference system	-
------------------	---

reference diaphragm	-
---------------------	---

reference electrolyte	-
-----------------------	---

bridge electrolyte	-
--------------------	---

electrode head	S7 (RedCap)
----------------	-------------

cable length	-
--------------	---

connector	-
-----------	---

cell constant (cm ⁻¹)	-
-----------------------------------	---

cell type	-
-----------	---

indicating metal	silver pin
------------------	------------

Storage	dry
----------------	-----

Special features

polypropylene shaft

silver pin 6 mm diameter, 10 mm length
--

Choice of cables

All Mettler-Toledo meters have BNC connectors. Recommended cable: BNC 1.0 m (59902392)

Our cable selector guide offers a large choice of S7 cables for other connectors. Red S7 connectors are compatible with grey S7 connectors. The color is only design-related.

Scale 1:1