

MS140-300

Portable Medical Electronic Handrail Scale



- Accuracy +/- 3 divisions
 (0.6 lb, 0.3 kg, 15 oz)
- Capacity
 – 660 lb x 0.2 lb; 300 kg x 0.1 kg; 660 lb x 5 oz
- Construction
 Stainless steel weighing surface with mild steel foam covered handrails. Rigid base design.
- Power–6 VDC, 500mA wall mount adapter (included)
- Weight Hold

 May be used to freeze and store the displayed weight so that you may attend to the patient's needs
- Integrated Wheel Kit- For easy portability



MS140-300 Portable Medical Electronic Handrail Scale

Features

Simplicity of design and rugged construction are the key features of the MS140-300 Portable Medical Scale.

Featuring a 520 mm (L) x 520 mm (W) / 20.5" x 20.5", low-profile (65 mm (H) / 2.5" H) platform with sturdy handrails and padded grips, the MS140-300 digital scale was designed to be the perfect choice for clinics, health clubs, or even home use.

This medical scale has the ability to print weight data offering a real time saving device for doctor's offices, long-term care facilities, and other medical institutions.

Heavy-duty, yet lightweight & portable this scale is suitable for almost any application requiring accuracy and medium to heavy weight capacity. The repeatability is unbeatable and will offer years of continued usage. It comes with an easy-to-read LCD backlit display and operates on AC power. This medical scale is durable and reliable at an affordable price.



Specifications

Capacity (lb) Capacity (kg) Capacity (lb/oz) 660 lb x 0.2 lb 300 kg x 0.1 kg 660 lb x 5 oz

Accuracy: +/- 3 divisions (0.6 lb, 0.3 kg, 15 oz)

Platform Dimension:

Scale: 520 mm (L) x 520 mm (W) x 1,135 mm (H)

20.5" L x 20.5" W x 44.25" H

Weighing Platform:

520 mm (L) x 520 mm (W) x 65 mm (H)

20.5" L x 20.5" W x 2.5" H

Display: 17 mm / 0.625V high LCD

Indicator: 165 mm (L) x 30 mm (W) x 81 mm (H)

6.5" L x 1.2 W x 3.2" H

Serial Port: One (1) full duplex RS232 with cable Function Keys: Hold/Print, Units, Tare, On/Off/Zero

Unit: lb, kg, lb/oz

Power: 6 VDC, 500mA adapter (included) Optional Power Adapter: UK 816965002887

EU 816965002870

Operating Temp: 5° C to 35° C / 41° F to 95° F

RS-232: Supplied complete with 10' cable (9 pin to 9 pin type connector, bi-directional full duplex)

Product Weight: 19 kg / 42 lb

Packaging Dimension:

1168 mm (L) x 635 mm (W) x 152 mm (H)

46" L x 25" W x 6" H

Shipping Weight: 25 kg / 55 lb

Packaging Contents: Manual, platform, 3 piece handles, power supply, RS-232 cable, indicator, 3 hex wrenches

and pull through string

Approvals: Full FCC compliance, UL/CUL power adapter

Product Code: 816965004829



a division of Avery Weigh-Tronix, LLC

Brecknell

1000 Armstrong Dr., Fairmont, MN 56031 Toll free: 800-637-0529, Tel: 507-238-8702 Fax: 507-238-8271

e-mail: sales@brecknellscales.com www.brecknellscales.com



Brecknell

Foundry Lane, Smethwick West Midlands B66 2LP UK Tel: +44 (0) 845 246 6717 Fax: +44 (0) 845 246 6718 Email: sales@brecknellscales.co.uk www.brecknellscales.co.uk





