



BARIATRIC PORTABLE WHEELCHAIR SCALE

Light in Weight... Heavy in Capacity

The BRW1000's integral handles make it easy to lift and lower into position. Built-in wheels permit the scale to be rolled anywhere it's needed.

DETECTO'S BRW1000 Portable Bariatric Wheelchair Scale



FEATURING THE POWERFUL 758C WEIGHT INDICATOR

- Pushbutton and Keypad Tare
- Body Mass Index Calculation
- Selectable Key Lockout
- ID Storage With and Without a Known Tare Weight
- EMR-Ready Serial
- Connectivity



Batterypowered, so you can take the scale to the patient.

Detecto's "Roll-A-Weigh" BRW1000

The suprisingly lightweight design and built-in handles and wheels make Detecto's BRW1000 "Roll-A-Weigh" Scale capable of being moved wherever you need it for bariatric weighing. Whether's it's multiple weighing locations at a single facility or a stationary location, the "Roll-A-Weigh" gives you portability, ease of use and accuracy.

Designed for bariatric weighing "on-the-move", this portable scale features a 758C digital weight indicator that can be mounted on the scale, or remotely on a wall, or desk for convenient viewing. Heavy-duty construction, large rubber deck safety plate platform with integral ramp, and the 1,000 lb capacity makes this the perfect scale for weighing patients in manual or power wheelchairs, dialysis patients, or walk-ons.

SPECIFICATIONS

Platform Width: Power Source: Capacity: Display: 40" (accomodates bariatric wheelchairs) 6 "C" size Alkaline or Nickel-Cadium cells (not included) or optional AC adapter. 1,000 lb x 0.5 lb / 400 kg x 0.2 kg 0.7" high-contrast LCD

Model BRW1000



Detecto Scale reserves the right to improve, enhance, or modify features and specifications without prior notice.



Sold By: