


PH: (417) 673-4631 FAX: (417) 673-5001 WWW.DETECTOSCALE.COM

DR400-7586

DIGITAL PORTABLE PHYSICIAN SCALE



DR400-758C SPECIFICATIONS

Capacity: 400 lb x 0.2 lb / 180 kg x 0.1 kg

Platform Size: 12" W x 12" D x 2.25" H / 305 mm W x 305

D mm x 57 mm

Display: 0.5" / 12.7 mm high, liquid crystal,

reflective type

Keypad: 19 key, color-coded, sealed membrane

Height Range: 3 to 7 ft / 90 to 215 cm

Power: 6 "C" size Alkaline or NIckel-Cadmium

cells (not included) or optional AC adapter

Shipping Weight: 13 lb / 6 kg

Detecto's reliable DR400-758C digital physician scale is lightweight, portable, and battery powered, making it perfect for mobile clinics, home care nurses, or wherever a portable scale with Body Mass Index is needed. The unit has a slip-resistant, textured mat to ensure patient safety when getting on and off the scale.

The remote digital indicator comes with a wall-mount bracket and easily mounts on a counter top or wall. The indicator features a large 0.5" / 12.7 mm high-contrast LCD, selectable units, Body Mass Index, pushbutton and keypad tare, and automatic display lock. The DR400-758C quickly calculates Body Mass Index by simply entering the patient's height while weighing.

Optional Wall-Mount Height Rod

Model 3PHTROD-WM

Measures 30 - 78 in and 76 - 198 cm

Optional Black Carrying Case

Model DR400C-CASE

Carrying case with handle protects the DR400-758C scale for easy transport.






SPECIFICATORS ONE

DR400C Patient Scale

Accurate and reliable, Detecto's economical DR400C platform scale is lightweight and portable, perfect for mobile clinics and home care nurses. The remote indicator features a large 1" display, units conversion, and tare. Battery powered so you can have a scale anywhere you need it.



SPECIFICATIONS

Capacity: 400 lb x .5 lb / 181 kg x .2 kg

Display: 4 digit, 7 segment, liquid crystal display, 1" / 25 mm high

Keypad: 4 keys: lb/kg, Tare/Hold, On, Off

Platform: 12" square, painted steel, non-skid surface

Units: Pounds, kilograms

Load Cells: 1

Power: 9 volts DC supplied by 6 "AA" alkaline batteries or

115 volt AC adapter (included)


OPTIONAL CARRYING CASE AVAILABLE

Cardinal Scale reserves the right to improve, modify, or enhance features and specifications without prior notice.


PH: 417-673-4631 • FAX: 417-673-5001 P.O. BOX 151, WEBB CITY, MO 64870 USA

WEB SITE: http://www.detectoscale.com





GP-400-758C Portable Scale



The perfect scale for those who need a lightweight, portable scale.

The GP-400-758C, with its detachable indicator and carrying handle, is very lightweight. The large, easy-to-clean platform incorporates a slip-proof mat which, when combined with its nonskid plastic feet, assures stability.

This large capacity, general purpose scale, with its exceptionally sensitive graduations, features automatic zero, auto power off, and motion detector. The indicator has a large, easy-to-read LCD display and operates on conventional "C" size batteries. Instructions are printed on the keypad, making it simple to use.

Specifications

Capacity: 400 lb x 0.2 lb (180 kg x 0.1 kg)
Display: .7" large, easy-to-read LCD

• **Power:** Six (6) "C" size Alkaline batteries

• Features: 758C indicator

Polyethylene feet Slip-proof surface

• Net Weight: 14 lb

Platform Size: 18" x 14
Functions: Automatic zero

Motion detector

Auto power off

(selectable 1 to 9 minutes) On/off, zero, units, hold/

release