GBC Bench Counting Scale








Speed, Performance, Value - Perfect Balance


- Available February 2011
- Date and time
- Selectable multi-language text for printing
- · Preset counting to a value
- Large backlit LCD display
- Overload protection
- Non-slip adjustable levelling feet
- · Large stainless steel platform
- Pillar mounted
- · Full colour coded keypad with numeric entry
- Precision load cell technology

- Simple user-friendly operation
- Splashproof to protect from accidental spills
- Easy access to the rechargeable battery
- RS-232 bi-directional interface
- External calibration
- Rechargeable battery up to 90 hours operation
- Auto sleep / power down function to save battery life
- Low battery indication
- AC Adapter
- Counting
- Parts counting









www.adamequipment.com



Technical Specifications

Model	GBC 35a	GBC 70a	GBC 130a
Maximum capacity	35lb/16kg	70lb/32kg	130lb/60kg
Readability	0.001lb/0.5g	0.002lb/1g	0.005lb/2g
Repeatability (S.D.)	0.002lb / 1g	0.002lb / 1g	0.01lb / 4g
Linearity (±)	0.002lb / 1g	0.002lb / 1g	0.01lb / 4g
Typical	2-3 Sec	-	2-3 Sec
stabilisation time			
Pan size	11.8×15.7" / 300×400mm		

General Specifications

Power supply	12vDC 800mA adapter / 6v 4ah rechargeable battery	
Calibration	Automatic External	
Display	Large LCD backlit display	
Overall dimensions	33.9×15.7×24.4" / 860×400×620mm	
(wxdxh)		
Dimensions notes	Shipping Weight: 16Kg Dims: 90 x 46 x 21cm	

Key Accessories

- · RS-232 to USB Interface Cable (must be ordered with #3074010266 |
- AdamDU Data Collection Program
- · In-use Cover
- · Adam Printer
- · Thermal Paper for Adam Printer (5 Pack)
- · RS-232 Cable







Adam Equipment - offices worldwide

www.adamequipment.com