

**MORE
SAVINGS**

ON APPLICATION-DRIVEN
CENTRIFUGE PACKAGES

Details inside!



**Thermo Scientific Heraeus
General Purpose Centrifuges**

more samples, more speed
your answers faster

Featuring breakthrough Thermo Scientific technologies -



Easy and secure push-button Auto-Lock® rotor exchange for application versatility and cleaning convenience




Certified ClickSeal® bucket sealing system provides glove-friendly, one-handed snap-on covers to replace screw caps and clips



Thermo Scientific Fiberlite carbon fiber rotors provide lightweight and durable performance, with high g-force (up to 25,000)

Thermo
SCIENTIFIC

Choose a Thermo Scientific Heraeus general purpose centrifuge package!



CELL CULTURE

Outstanding throughput in compact footprint



BLOOD PROCESSING

Exceptional capacity with certified biocontainment

Heraeus® Multifuge® X1 and Megafuge® 16

16x50 ml tubes up to 5000 rpm/4696 xg	Cat No.	Price
Heraeus Megafuge 16 Centrifuge, 120 V	75318382	\$5304
Heraeus Megafuge 16R Centrifuge, 120 V	75718382	\$6985
Heraeus Multifuge X1 Centrifuge, 120 V	75118382	\$5967
Heraeus Multifuge X1R Centrifuge, 120 V	75518382	\$7291

PACKAGE INCLUDES

TX-400 Swinging Bucket Rotor
400 ml Round Buckets (set of 4)
Adapters 4 x 50 ml conical (set of 4)
Adapters 9 x 15 ml conical (set of 4)

Save up to
\$2500



36x5 ml tubes up to 5500 rpm/5500 xg	Cat No.	Price
Heraeus Megafuge 16 Centrifuge, 120 V	75310885	\$4233
Heraeus Megafuge 16R Centrifuge, 120 V	75710885	\$5915

PACKAGE INCLUDES

TX-200 Swinging Bucket Rotor
200 ml Round Buckets (set of 4)
ClickSeal Biocontainment Lids for Round Buckets (set of 4)
Adapters 7 x 15/10 ml blood collection (set of 4)
Adapters 8 x 5/7 ml blood collection (set of 4)

Save up to
\$2200

76x5 ml tubes up to 5000 rpm/4700 xg	Cat No.	Price
Heraeus Megafuge 16 Centrifuge, 120 V	75318180	\$5304
Heraeus Megafuge 16R Centrifuge, 120 V	75718180	\$6985

PACKAGE INCLUDES

TX-400 Swinging Bucket Rotor
400 ml Round Buckets (set of 4)
ClickSeal Biocontainment Lids for Round Buckets (set of 4)
Adapters 14 x 15/10 ml blood collection (set of 4)
Adapters 19 x 5/7 ml blood collection (set of 4)

Save up to
\$2700

Heraeus Multifuge X3 and Megafuge 40

28x50 ml tubes up to 4700 rpm/4800 xg	Cat No.	Price
Heraeus Megafuge 40 Centrifuge, 120 V	75043839	\$7597
Heraeus Megafuge 40R Centrifuge, 120 V	75193839	\$9636
Heraeus Multifuge X3 Centrifuge, 120 V	75013839	\$8464
Heraeus Multifuge X3R Centrifuge, 120 V	75163839	\$10146
Heraeus Multifuge X3F Centrifuge, 120 V	75313839	\$9475
Heraeus Multifuge X3FR Centrifuge, 120 V	75373839	\$11156

PACKAGE INCLUDES

TX-750 Swinging Bucket Rotor
750 ml Round Buckets (set of 4)
Adapters 7 x 50 ml conical (set of 4)
Adapters 14 x 15 ml conical (set of 4)

Save up to
\$4000

40x50 ml tubes up to 3500 rpm/2700 xg	Cat No.	Price
Heraeus Megafuge 40 Centrifuge, 120 V	75047473	\$7036
Heraeus Megafuge 40R Centrifuge, 120 V	75197473	\$9024

PACKAGE INCLUDES

BIOLiner® Swinging Bucket Rotor
BIOLiner Buckets (set of 4)
Adapters 10 x 50 ml conical (set of 4)
Adapters 22 x 15 ml conical (set of 4)

Save up to
\$3800

108x5 ml tubes up to 4700 rpm/4800 xg	Cat No.	Price
Heraeus Megafuge 40 Centrifuge, 120 V	75041923	\$7597
Heraeus Megafuge 40R Centrifuge, 120 V	75191923	\$9636

PACKAGE INCLUDES

TX-750 Swinging Bucket Rotor
750 ml Round Buckets (set of 4)
ClickSeal Biocontainment Lids for Round Buckets (set of 4)
Adapters 21 x 10 ml blood collection (set of 4)
Adapters 27 x 5/7 ml blood collection (set of 4)

Save up to
\$4000

196x5 ml tubes up to 3500 rpm/2700 xg	Cat No.	Price
Heraeus Megafuge 40 Centrifuge, 120 V	75047271	\$7087
Heraeus Megafuge 40R Centrifuge, 120 V	75197271	\$9126
Heraeus Multifuge X3F Centrifuge, 120 V	75317271	\$9475
Heraeus Multifuge X3FR Centrifuge, 120 V	75377271	\$11717

PACKAGE INCLUDES


BIOLiner Swinging Bucket Rotor
BIOLiner Buckets (set of 4)
ClickSeal Biocontainment Lids for BIOLiner Buckets (set of 4)
Adapters 37 x 10 ml blood collection (set of 4)
Adapters 49 x 5/7 ml blood collection (set of 4)

Save up to
\$3800



Biocontainment certification by HPA, Porton Down, UK (formerly CAMR).

Designed around you and your applications!



MICROPLATE PROCESSING

Incredible productivity in sealed conditions




CARBON FIBER

Super speed performance on the bench

6 standard or 2 deepwell up to 4000 rpm/2272 xg	Cat No.	Price
Heraeus Megafuge 16 Centrifuge, 120 V	75312452	\$4029
Heraeus Megafuge 16R Centrifuge, 120 V	75712452	\$5915



PACKAGE INCLUDES

 M-20 Microplate Rotor

Save up to
\$3200

8x50 ml up to 24,400 xg	Cat No.	Price
Heraeus Multifuge X1R Centrifuge, 120 V	75512963	\$8617

PACKAGE INCLUDES

 TX-400 Swinging Bucket Rotor
400 ml Round Buckets (set of 4)
 Fiberlite® F15-8x50c Rotor

Save up to
\$3300



16 standard or 4 deepwell up to 4700 rpm/4800 xg	Cat No.	Price
Heraeus Multifuge X3 Centrifuge, 120 V	75010817	\$7903
Heraeus Multifuge X3R Centrifuge, 120 V	75190817	\$9781

PACKAGE INCLUDES



TX-750 Swinging Bucket Rotor
Carriers for Microplates
and T-75 Horizontal Flasks

Save up to
\$4600



6x250 ml up to 18,500 xg	Cat No.	Price
Heraeus Multifuge X3R Centrifuge, 120 V	75160062	\$10953

PACKAGE INCLUDES



 TX-750 Swinging Bucket Rotor
750 ml Round Buckets (set of 4)
 Fiberlite F14S-6x250LE Rotor

Save up to
\$5000

8x50 ml up to 24,400 xg

Heraeus Multifuge X3R Centrifuge, 120 V	75160063	\$10961
---	----------	---------

PACKAGE INCLUDES

 TX-750 Swinging Bucket Rotor
750 ml Round Buckets (set of 4)
 Fiberlite F15-8x50c Rotor

Save up to
\$5000



and maintain the
flexibility of your
swing-out

* 50 ml (75003683) and 15 ml (75003682) conical adapters for TX-400 rotor sold separately.
** 50 ml (75003685) and 15 ml (75003684) conical adapters for TX-750 rotor sold separately.



Your centrifuge is essential to your work. It can also be a bottleneck to efficient research. Unless you choose the family of Thermo Scientific Heraeus general purpose centrifuges designed around you and your applications to deliver flexible, efficient sample processing.



Heraeus® Multifuge® X1 and Multifuge X3 centrifuge series

With the latest user interface and extensive rotor selection providing 23+ rotor configurations, meet all processing needs at RCF up to 25,000 xg.

Heraeus Megafuge® 16 and Megafuge 40 centrifuge series

Accelerate your research by processing 4 to 40 x 50 ml conical tubes or from 32 to 196 blood tubes with our extensive selection of biocontainment rotor solutions, all requiring minimum bench space.

Also available in compact floor-standing configuration – Heraeus Multifuge X3F centrifuge series.

**ACT
NOW!**

**For more details,
contact your sales representative
or call 866 984 3766.**

Special pricing on Heraeus Multifuge X1, Heraeus Multifuge X3, Heraeus Multifuge X3F, Heraeus Megafuge 16 and Heraeus Megafuge 40 packages expires on December 31, 2012. Valid in North America only. Void where prohibited. Not valid with any other offer, contract or agreement.

©2012 Thermo Fisher Scientific Inc. All rights reserved. HERAEUS is a registered trademark of Heraeus Holding GmbH licensed to Thermo Fisher Scientific. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries.

North America: USA/Canada +1 866 984 3766 (866-9-THERMO)

Europe: Austria +43 1 801 40 0, Belgium +32 2 482 30 30, France +33 2 2803 2180, Germany national toll free 08001-536 376, Germany international +49 6184 90 6940, Italy +39 02 02 95059 434-254-375, Netherlands +31 76 571 4440, Nordic/Baltic countries +358 9 329 100, Russia/CIS +7 (812) 703 42 15, Spain/Portugal +34 93 223 09 18, Switzerland +41 44 454 12 12, UK/Ireland +44 870 609 9203

Asia: China +86 21 6865 4588 or +86 10 8419 3588, India toll free 1800 22 8374, India +91 22 6716 2200, Japan +81 45 453 9220,

Other Asian countries +852 2885 4613 Countries not listed: +49 6184 90 6940 or +33 2 2803 2180

www.thermoscientific.com/centrifuge

Thermo
SCIENTIFIC

Part of Thermo Fisher Scientific