

Capacity: 400 lb



### Advantages

- ✓ **400 lb Capacity Floor scale**  
The high capacity of this floor scale provides more versatile use
- ✓ **8" Raised Dial Display**  
Large, easy-to-read dial display is raised above the scale base to allow for a larger platform
- ✓ **Weighs in LB**  
Dial display shows weight in pounds

**Health o meter®**  
P r o f e s s i o n a l

Capacity: 400 lb



Model	Master Pack	Shipping Dims. (L x W x H) / Cube	Weight (LB)	UPC Code
160LB	2	17.25" x 15.25" x 22.13" / 3.37	22 lb	892076002053

**Specifications**

- Capacity: 400 lb
- Graduation: 1 lb
- 8" dial (raised 4")
- Platform Size: 11-1/2" (w) x 13" (d)
- 1 Year Limited Warranty