



Advantages

- ☒ **EMR Connectivity**
 Scale has the ability to interface with a computer, EMR software or other electronic device via various Health o meter® Professional Connectivity Solutions
- ☒ **Calculates Body Mass Index (BMI)**
 Scale functions include the ability to calculate the patient's Body Mass Index (BMI)
- ☒ **Value-Added Features**
 Operating functions (LB/KG Conversion, LB/KG Lock Out, BMI, Zero Out/Tare, Hold/Release, Reweigh, Recall, Auto Zero, Auto Off) included at no additional cost
- ☒ **Hospital Grade Construction**
 Durable steel structure with large sturdy stretcher enables patients to be weighed accurately and comfortably



Model	Master Pack	Shipping Dims. (L x W x H) / Cube (3 cartons)	Weight (LB)	UPC Code
2000KL	Bar	68.375" x 3.25" x 2.125" / 0.27	9.25 lb	892076002466
	Hoyer Lift	51.625" x 25.875" x 13.125" / 10.15	76 lb	
	Assembly	32.375" x 20.375" x 6.375" / 2.43	45 lb	



Specifications

- Capacity: 400 lb / 181 kg
- Resolution: 0.2 lb / 0.1 kg
- Connectivity via USB
- 3/4" LCD display
- Functions: LB / KG Conversion, LB / KG lock Out, Body Mass Index (BMI), Zero Out / Tare, Hold / Release, Reweigh, Auto Zero, Auto Off
- Canvas Size: 70½" (w) x 30" (d) / 1791 mm (w) / 762 mm (d)
- (4) Wheels for ease of movement, 1 locking wheel for safety
- Power Source: 110-240V adapter included, 6 AA batteries (not included)
- 2 Year Limited Warranty