

**DISPLAYS WEIGHT IN POUNDS, KILOGRAMS & GRAMS**

**Health o meter®**  
P r o f e s s i o n a l



## Advantages

- Motion Detection**  
Motion sensing technology compensates for baby's movement and high-resolution readings provide more accurate weight measurements, ideal for the neonatal environment
- Convenient Features**  
Plastic tray is easy-to-clean and heavy-duty built-in handles allow for portability
- Built-In Measuring Tape**  
Easy-to-read measuring tape allows user to measure both weight and height simultaneously
- Industry Leading 3 Year Warranty**  
Health o meter® Professional's industry leading warranty protects you against defects of materials or workmanship
- Value-Added Features**  
Operating functions (LB/KG/GM Lockout, LB/GM Conversion, LB/KG Conversion, Zero Out/Tare, Reweigh/Recall, Auto Zero, Auto Off) included at no additional cost

**Health o meter®**  
P r o f e s s i o n a l  
[www.hom scales.com](http://www.hom scales.com)

**HIGH RESOLUTION  
PORTABLE SCALE**



**NON-CONDUCTIVE  
WHEELS**

**2210CART**

Model	Shipping Dims. (L x W x H) / Cube	Weight (LB)	UPC Code
2210KL	27.875" x 18.625" x 7.625" / 2.29	30.5 lb	892076002329
2210CART	2 Cartons: 26.5" x 18.5" x 9" / 2.55	23 lb	
	22" x 17.75" x 2.5" / 0.56	17 lb	

## 2210KL



### 2210KL Specifications

- Capacity: 45 lb / 20 kg / 20,000 gm
- Resolution: 0 – 11300 g (0 – 25 lb) (0 - 10 kg) +/- 1g (+/- 0.1 oz)  
> 11300 g (> 25 lb) (> 10 kg) +/- 5 g (+/- 0.1 oz)
- 1 gram accuracy standard
- 7/8" LCD display
- Functions: LB / KG / GM Lockout, LB / GM Conversion, LB / KG Conversion, Zero Out / Tare, Reweigh / Recall, Auto Zero, Auto Off, Motion Detection
- Removable curved clear acrylic tray
- Tray Size: 24" (w) x 15-3/8" (d) x 3" (h) / 610 mm (w) x 391 mm (d) x 76 mm (h)
- Measuring Tape: 0" - 23" / 0 cm - 60 cm
- Heavy duty handles for ease of moving
- 6 C-cell batteries included
- 3 Year Limited Warranty

### 2210CART Specifications

- Capacity: 200 lb / 90 kg
- Top shelf has handle for easy grip
- Non-conductive wheels / 2 wheels with locks
- Cart Dimensions: 26-1/2" (w) x 18" (d) x 34-1/8" (h) / 673 mm (w) x 457 mm (d) x 867 mm (h)
- Lockable Cabinet: 23" (w) x 16-3/4" (d) x 17" (h) / 584 mm (w) x 427 mm (d) x 432 mm (h)
- Storage Shelf is easy to clean, durable blue plastic
- Use with pediatric scales: 2210KL, 553KL, 549KL
- 1 Year Limited Warranty