


conductivity



Scale 1:1

Order number 51344110

former order number -

Specifications

pH range -

zero point (pH) -

temperature range 0 - 100 °C

temperature probe NTC 30 kΩ

shaft material epoxy

membrane glass -

membrane resistance -

reference system -

reference diaphragm -

reference electrolyte -

bridge electrolyte -

electrode head -

cable length 1.8 m

connector Mini-DIN (IP67)

cell constant (cm⁻¹) 0.57

cell type 4 graphite poles

indicating metal -

Storage dry**Additional specifications**

measuring range 0.01 - 1000 mS/cm

Intelligent Sensor Management (ISM)

IP67

Instrument compatibility

SevenGo Duo™ SG23 / SG78