


conductivity



Order number 51344126

former order number

Specifications

pH range -

zero point (pH) -

temperature range 0 - 100 °C

temperature probe NTC 30 kΩ

shaft material stainless steel

membrane glass -

membrane resistance -

reference system -

reference diaphragm -

reference electrolyte -

bridge electrolyte -

electrode head -

cable length 1.8 m

connector Mini-DIN (IP67)

cell constant (cm⁻¹) 0.105

cell type 2 steel poles

indicating metal -

Storage dry**Additional specifications**

measuring range 0.001 - 500 µS/cm

IP67

Instrument compatibility

SevenGo™ SG3 / SG7 and

SevenGo Duo™ SG23 / SG78

replaces the

InLab®740 SG

Scale 1:1