



Order number 51343195

former order number -

Specifications

pH range 0 - 14

zero point (pH) 7.00 ± 0.25

temperature range 0 - 100 °C

temperature probe -

shaft material glass

membrane glass HA

membrane resistance $< 600 \text{ M}\Omega$

reference system -

reference diaphragm -

reference electrolyte -

bridge electrolyte -

electrode head S7 (grey)

cable length -

connector -

cell constant (cm^{-1}) -

cell type -

indicating metal -

Storage 3 mol/L KCl

Choice of cables

All Mettler-Toledo meters have BNC connectors.
Recommended cable: BNC 1.2m (52300004)

Our cable selector guide offers a large choice of S7 cables for other connectors. Grey S7 connectors are compatible with red S7 connectors. The color is only design-related.

replaces the

InLab®201