



**Order number** 51343196

former order number -

#### Specifications

pH range 0 - 12

zero point (pH)  $7.00 \pm 0.25$

temperature range 0 - 130 °C

temperature probe -

shaft material glass

membrane glass A41, thick-walled

membrane resistance  $< 700 \text{ M}\Omega$

reference system -

reference diaphragm -

reference electrolyte -

bridge electrolyte -

electrode head S7 (grey)

cable length -

connector -

cell constant ( $\text{cm}^{-1}$ ) -

cell type -

indicating metal -

**Storage** 3 mol/L KCl

#### Special features

extra strong glass membrane

maximum temperature 130 °C

#### Choice of cables

All Mettler-Toledo meters have BNC connectors.  
Recommended cable: BNC 1.2m (52300004)

Our cable selector guide offers a large choice of S7 cables for other connectors. Grey S7 connectors are compatible with red S7 connectors. The color is only design-related.

replaces the

**265-VM-S7/120**