


ORP / redox



Scale 1:1

Order number 51343205

former order number -

Specifications

pH range -

zero point (pH) -

temperature range 0 - 100 °C

temperature probe -

shaft material glass

membrane glass -

membrane resistance -

reference system ARGENTHAL™, Ag⁺ trap

reference diaphragm ceramic

reference electrolyte 1 mol/L KNO₃

bridge electrolyte -

electrode head S7 (grey)

cable length -

connector -

cell constant (cm⁻¹) -

cell type -

indicating metal silver ring

Storage 1 mol/L KNO₃**Special features**

silver ring

standard electrolyte 1 mol/L KNO₃ARGENTHAL™ reference system with Ag⁺ trap**Choice of cables**

All Mettler-Toledo meters have BNC connectors.
Recommended cable: BNC 1.2m (52300004)

Our cable selector guide offers a large choice of S7 cables for other connectors. Grey S7 connectors are compatible with red S7 connectors. The color is only design-related.

replaces the

InLab®503