



**Order number** 51343192

former order number -

#### Specifications

pH range -

zero point (pH) -

temperature range 0 - 130 °C

temperature probe -

shaft material glass

membrane glass -

membrane resistance -

reference system ARGENTHAL™

reference diaphragm triple ceramic

reference electrolyte 3 mol/L KCl

bridge electrolyte -

electrode head S7 (grey)

cable length -

connector -

cell constant (cm<sup>-1</sup>) -

cell type -

indicating metal -

**Storage** 3 mol/L KCl

#### Special features

triple ceramic diaphragm

electrolyte refill hose coupling

#### Choice of cables

All Mettler-Toledo meters have BNC connectors.  
Recommended cable: BNC 1.2m (52300004)

Our cable selector guide offers a large choice of S7 cables for other connectors. Grey S7 connectors are compatible with red S7 connectors. The color is only design-related.

replaces the

**360-T-S7/120**