



Order number 51343054

former order number -

Specifications

pH range 0 - 14

zero point (pH) 7.00 ± 0.25

temperature range 0 - 100 °C

temperature probe NTC 30 kΩ

shaft material glass

membrane glass HA

membrane resistance < 600 MΩ

reference system ARGENTHAL™, Ag⁺ trap

reference diaphragm ceramic

reference electrolyte 3 mol/L KCl

bridge electrolyte -

electrode head MultiPin™

cable length -

connector -

cell constant (cm⁻¹) -

cell type -

indicating metal -

Storage 3 mol/L KCl

Special features

integrated NTC 30 kΩ temperature probe

ARGENTHAL™ reference system with Ag⁺ trap

Choice of cables

All Mettler-Toledo meters have BNC connectors.
Recommended cable: BNC/RCA (Cinch) 1.2 m
(52300009)

Our cable selector guide offers a large choice
of MultiPin™ cables for other connectors.

replaces the

InLab®410