

# **PGW Precision Balances**



# "The Affordable Choice"

www.adamequipment.com



# **PGW Precision Balances**

## "The Affordable Choice"

The PGW Series of Precision Balances are designed to meet the needs of Laboratories, Schools, Industry and commercial users offering a wide range of applications.

- RS-232 bi-directional interface
- Second Communication port
- Kensington<sup>™</sup> type lock bracket
- External calibration
- Internal calibration option
- GLP print outs
- Date and time
- Colour coded keys

- Backlit LCD with dual display
- Up to 18 weighing units (including one custom unit)
- Windshield supplied as standard with models with 0.001g readability
- Capacity tracker
- Multilingual display
- PGW PGW PGW PGW Model PGW PGW PGW PGW 153e 253e 453e <u>7</u>53e 1502e 2502e 3502e 4502e Maximum 150g 250g 450g 750g 1500g 2500g 3500g 4500g capacity Readability 0.001g 0.01g Full Tare range Repeatability 0.001g 0.01g Linearity (±) 0.002g 0.02g Units of measure gram, milligram, kilogram, carat, grain, newton, dram, pound, ounce, troy ounce, pennyweight, momme, tael (Taiwan, Hong Kong and Singapore), tola, tical and 1 custom unit Bi-directional RS-232 (Standard) 15 VDC, 50/60 Hz, 800 mA Interface Power supply Internal or External, Selectable automatic calibration due to change Calibration in time or temperature Backlit LCD with dual digits (24 mm high) and capacity tracker Display Draught shield Supplied as standard on models with 1mg readability Die cast aluminium housing (glass weighing chamber for models Housing with 1mg readability) 7.6×7.6" / 192×192mm 5.5×5.5" / 140×140mm Pan size 9.9×14.1×4.1" / 251×358×104mm Overall dimensions (wxdxh) Applications Weighing, Net / total, Check weighing, Percentage weighing, Weight accumulation, Averaging, Parts counting, Animal / dynamic weighing, Density determination, Below balance weighing





#### **Optional Accessories:**

- Below Balance Hanger
- Density Determination Kit
- (mg models only)
- In-use Cover (mg models only)
- Dust Cover
- Security lock and cable
  In-use Cover



### Adam Equipment have offices in the UK, USA, South Africa and China

## www.adamequipment.com

SPECIFICATIONS