

STAINLESS STEEL

Available in the world's standard, the Eye-Level Mechanical Beam Physician Scale and also in the state-of-the-art, EMR-ready Waist-High Digital Scale.



UNIQUELY DETECTO - UNIQUELY STAINLESS STEEL



Model 339S

Model 8437S



• **Hygiene** The time-proven benefits of *stainless*

steel - cleans easily and completely.

• **Beauty** Shiny, attractive **stainless steel** - fits

well in all modern decors.

• Protection Rust-resistant *stainless steel* protects

the scale - especially useful in high

humidity areas.



Stainless Steel MECHANICAL Health Care Scales

Enjoy the classic elegance of stainless steel in your Detecto mechanical weighbeam scale. A platform cover that lifts off easily combined with the stainless steel casing makes cleanup a breeze. Water damage from mopping and high humidity areas around swimming pools is a thing of the past with Detecto's stainless steel health care scales. Scale height is 59" (1.5 m) and the die-cast beam can be easily read from either side.

Model	Capacity		Height Rod	
	400 lb x 4 oz	180 kg x 100 g	30 - 78 in 76 - 198 cm	Shipping Weight
439S	•		•	38 lb / 17.2 kg
2391S		•	•	37 lb / 16.7 kg
437S	•			35 lb / 15.8 kg
2371S		•		35 lb / 15.8 kg
339S	•	•	•	38 lb / 17.2 kg
337S	•	•		35 lb / 15.8 kg

Now, for the first time, you can combine the classic Detecto physician scale with the time-honored beauty and protection of stainless steel.

Stainless Steel DIGITAL Health Care Scales

With an ultra-modern design, Detecto's stainless steel digital health care scales offer an extra-tough finish which stands up to water damage from mopping around the base and high-humidity areas such as locker rooms. Enjoy the power of digital weighing with the convenience of "at-a-glance" weight readouts and an AC or battery-powered LCD display. With the push of a button, readings appear in pounds OR kilograms.

Model	Capacity 400 lb x 0.2 lb / 180 kg x 0.1 kg	Height Rod 30 - 78 in 76 - 198 cm	Shipping Weight
8430S	•	•	34 lb / 15.4 kg
8437S	•		31 lb / 14 kg



 $Detecto Scale \ reserves \ the \ right to \ improve, enhance \ or \ modify \ features \ and \ specifications \ without \ prior \ notice. All \ registered \ trademarks \ are \ the \ property \ of \ their \ respective \ owners.$



