


**ThermoFisher**  
SCIENTIFIC  
275 AIKEN ROAD  
ASHEVILLE, NC 28804

THE INFORMATION CONTAINED HEREIN IS  
THE PROPRIETARY DATA OF THERMO  
FISHER SCIENTIFIC. IT IS NOT TO BE  
REPRODUCED OR DISCLOSED IN WHOLE  
OR IN PART WITHOUT PRIOR WRITTEN  
CONSENT OF THERMO FISHER SCIENTIFIC.  
DO NOT SCALE THIS DRAWING.

TITLE **REFRIGERATOR  
GLASS DOOR,45 CuFt**

-	60088	2011-04-05	INITIAL RELEASE
REV.	ECN #	DATE	DESCRIPTION
REVISIONS			

3rd ANGLE  
PROJECTION



DWG. NO.


SIZE **S**

SHEET 1 OF 2

**MR45SS-GAEE**

## GENERAL SPECIFICATIONS

1. EXTERIOR FINISH: STAINLESS STEEL FRONT, PAINTED GRAY SIDES
2. EXTERIOR DIMENSION:  
WIDTH 1338.6mm(52.7") DEPTH 812.8mm(32") HEIGHT 2057.1mm(81")
3. INTERIOR FINISH: STAINLESS STEEL
4. INTERIOR DIMENSION:  
WIDTH 1244.6mm(49") DEPTH 685.8mm(27") HEIGHT 1568.5mm(61.75")
5. INSULATION: FOAMED-IN-PLACE, CFC AND HCFC-FREE POLYURETHANE
6. DOOR: DOUBLE PANE, TEMPERED GLASS, SLIDING GLASS DOUBLE DOOR, KEY LOCK
7. SHELVES: GRAY-8 WIRE, EPOXY COATED, ADJUSTABLE IN 1/2" INCREMENTS
8. CASTERS: 4 STEM MOUNT
9. CONTROLLER: DIGITAL, FACTORY SET TO 4°C (39°F)
10. DEFROST: N/A
11. COMPRESSOR: 1/3 HP, REFRIGERANT R 134A
12. ELECTRICAL DATA: 115V 60Hz / 8.8 AMPS
13. PLUG TYPE: NEMA 5-15P
14. TOTAL CAPACITY: 45 Cu.ft
15. TEMPERATURE RANGE: 1 TO 12°C
16. AVERAGE TEMPERATURE STABILITY:  $\pm 3^{\circ}\text{C}$
17. APPROVALS: UL/cUL
18. SHIPPING WEIGHT: 255KG (561 LBS).

				<b>ThermoFisher</b> <b>SCIENTIFIC</b> 275 AIKEN ROAD ASHEVILLE, NC 28804		THE INFORMATION CONTAINED HEREIN IS THE PROPRIETARY DATA OF THERMO FISHER SCIENTIFIC. IT IS NOT TO BE REPRODUCED OR DISCLOSED IN WHOLE OR IN PART WITHOUT PRIOR WRITTEN CONSENT OF THERMO FISHER SCIENTIFIC. DO NOT SCALE THIS DRAWING.
				TITLE <b>REFRIGERATOR GLASS DOOR, 45 CuFt</b>		
-	-	-	SEE SHEET 1	3rd ANGLE PROJECTION 	DWG. NO.	
REV.	ECN#	DATE	DESCRIPTION		MR45SS-GAEE	
REVISIONS				SIZE <b>S</b>	SHEET 2 OF 2	