

# Precision Health Scale









# Precision Health Scale for Home, Athletic and Medical Use

# 10 Times More Accurate than Conventional Bathroom Scales

(350lb x 0.1lb)

The UC-321S Precision Scale is a sleek, stylish unit that offers more precise readings and functionality than those offered from traditional scales. This scale provides your weight to a tenth of a pound for precise information.

#### **Precise Measurements**

The UC-321S Precision Scale provides ultra precise weight results. Even the slightest increase or decrease in weight from an object like a glass of water registers on this scale. The UC-321S scale is extremely useful for diabetics and athletes who need to closely monitor and control their body weight within a certain range.



#### Large Easy-to Read Display

Superior Accuracy

Accurate to 0.1 lb up to 350 lbs.

#### **Target Weight Function**

Automatically calculates difference between desired target and actual weight

#### 31 Measurement Memory

Calculates difference between previous and current weight

#### Body Mass Index (BMI)

Automatically calculates a BMI value

#### Simple One Button Control

Easy toe operation, no need to bend over.

Specifications		
Capacity	350lb	
Resolution	0.1lb	
Display	LCD, character height, 1" height (25mm)	
Power	4 type AA batteries	
Battery Life	Approx. 2,000 measurements	
Weight	Approx. 4.85lb (2.2kg)	
Dimensions	12.6 x 12.3 x 1.4" (W x D x H)	
Auto Power Off	5 seconds	
Operating Temperature	50~95°F (10~35°C)	
RS-232C Computer Interface	Model UC-321PS only	
Warranty	3 years	



1555 McCandless Drive, Milpitas, CA 95035 (800) 726-3364 (408) 263-5333 Fax: 408-263-0119 Fax-On-Demand: 800-726-7099, press 2

Email: scales@andweighing.com www.andweighing.com





# Thin and Light

The UC-321S contains A&D's state-of-the-art technology in a stylish and low profile body. The UC-321S is one of the thinnest and lightest personal scales you can buy. Measuring only 1.4" thick, this scale is easy to step on and store.

## Target Weight Setting & 31 Built-in Memory Functions

With the UC-321S, you can set a target weight and it will tell you the difference between your actual and target weights. Additionally, you can store up to 31 readings in memory and recall the stored values by pressing a button.

### Body Mass Index Calculation

Body mass Index (BMI) takes into consideration an adult's weight-to-height ratio to determine if a person is overweight. The BMI value is calculated automatically and shown on the display.

The UC-321PS model includes the RS-232C computer interface for data collection into Windows™ software applications.